Instructions for Use

Page 2

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Seite 12

Mode d'emploi

Page 22

Istruzioni d'uso


Pagina 32

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Página 42

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Página 52

Instructions for Use

Introduction

Please read these Instructions for Use carefully and familiarize yourself with the application techniques.

The prefabricated SR Vivo TAC / SR Ortho TAC tooth line is a diagnostic, radiopaque tool that is designed to help dental technicians fabricate implant-supported restorations. Because of their special composition, these teeth are radiopaque.

Advantages

- A defined radiopacity achieved through industrial fabrication
- Moulds corresponding to those of the popular SR Vivodent / SR Orthotyp and SR Vivodent PE / SR Orthotyp PE tooth lines
- Easily prepared diagnostic aid
- The shape of radiopaque teeth can be modified with radiopaque monomers and polymers
- Strong bond with heat-curing as well as cold-curing polymers
- Prefabricated radiopaque teeth save time in the fabrication of stents
- The teeth ensure the accurate placement of dental implants during the planning phase
- Homogeneous dispersion of the radiopaque material in the teeth

Composition

The resin teeth are composed of the following ingredients:

Polymethyl methacrylate and barium sulphate
66–67 wt%

Methyl methacrylate and dimethacrylate 33–34 wt%

Additional components: catalysts and stabilizers
0.5 wt%

Modifier Monomer:

Methyl methacrylate and dimethacrylate 99 wt%
Catalysts and stabilizers 1 wt%

Modifier Polymer:

Polymethyl methacrylate and barium sulphate

> 99 wt%

Catalyst < 1 wt%

Indication

For the fabrication of fixed and removable implant-supported dental restorations in situations where the set-up of teeth is beneficial to the pre-surgical radiological diagnosis (CT, telerradiography as well as OPG), eg:

- severe dysgnathia
- severe atrophy of the alveolar ridge
- reconstructive surgery of the mouth, jaw or face

Pre-surgical planning of the restoration enables the ideal placement, location and angulation of the implant to be established.

Contraindication

The product is contraindicated for all the applications that are not specifically recommend by Ivoclar Vivadent as well as for permanent restorations in the mouth of patients.

In case of fixed metal prosthetics, metal fillings or metal osteosynthetic materials, the image quality of the CT and thus the recognizability of the radiopaque tooth may be impaired.

Side effects

To date, no systemic side effects are known. In rare cases, allergic reactions have been reported in conjunction with methyl methacrylate and polymethyl methacrylate materials. Do not use SR Vivo TAC / SR Ortho TAC or Modifier Polymer and Monomer if a patient is known to be allergic to any of the products' ingredients.

Processing

The radiopaque teeth SR Vivo TAC / SR Ortho TAC are set up according to conventional principles. The moulds of the radiopaque teeth are identical to those of the SR Vivodent / SR Orthotyp tooth line. SR Vivo TAC / SR Ortho TAC exhibit excellent, consistent radiopacity. Consequently, the two and three-dimensional relationship between the position of the implant and that of the subsequent superstructure can be established.

Fabrication of diagnostic stents

1. Laboratory

Diagnostic stents are fabricated to allow the implant-supported complete denture to be set up according to the instructions of the dentist, taking aesthetic, phonetic and functional aspects into account. In the process, SR Vivo TAC / SR Ortho TAC are used like conventional denture teeth.


The use of set-up aids such as the 2D/3D setting up template (Ivoclar Vivadent AG) is recommended.

Suggested procedure

Once the radiopaque teeth have been set up in wax, excess wax in the vestibular region is removed to lengthen the cervical region with the Monomer/Polymer Modifier until the base of the tooth rests on the model. Consequently, the entire length of the tooth is made of radiopaque material.

Wax model with lengthened cervical regions

Completed wax model

Completed model with lengthened cervical regions

Adjustment of shape

- Following the set-up and wax-up procedures, a key is fabricated using silicone or plaster.

Key of teeth set up in wax

- Remove wax from the model and teeth without leaving any residue. If a plaster key is used, the model and the key have to be isolated with Ivoclar Vivadent Separating Fluid. In addition, the tooth surfaces requiring adjustments must be roughened.

Reposition the teeth in the key

Check the accuracy of fit

- Mix the Polymer and Monomer Modifier in a 1.3 : 1 ratio and leave to rest uncovered for 2–3 minutes.

Mix the Monomer and Polymer Modifier

Pour the liquid material into the key

Secure the key

- Polymerization is conducted in a pressure device or in a heat/pressure polymerization unit (eg Ivomat), 15 minutes at 40–50 °C / 104–122 °F and 2–6 bar pressure.

After polymerization

If necessary, radiographic balls or similar aids may be integrated into the diagnostic stents where they are needed.

Finishing

Conventional finishing and polishing instruments are used for finishing (according to the conventional PMMA working techniques).

After finishing

Stent design for gaps


If requested by the dentist, a radiopaque alveolar ridge with mucosa can be created by coating this area with the Monomer/Polymer Modifier.

Fabrication of base plates

In order to reinforce or stabilize the diagnostic stents, a base plate is fabricated for the palatal or lingual region. It is made of clear or pink denture base material (eg ProBase Cold pink/clear) or a heat-curing polymer (eg ProBase Hot pink/clear).

ProBase Cold clear (also available in pink)

Space retainer for radiopaque Modifier material

When fixed restorations are planned (eg with an interdental gap), the base plate rests on the mesial or distal remaining dentition to stabilize and position the restoration.

Polymerized stents

Polishing

As the stents are used as a diagnostic tool and/or an aid during the surgical intervention, thorough prepolishing with goat's hair brushes and pumice as well as high gloss polishing with cotton buffers and polishing paste (Universal Polishing Paste from Ivoclar Vivadent) is recommended.

Check the occlusion

Finished diagnostic stent

Option

Diagnostic stent for complete maxillary denture

Transformation of the diagnostic stent into a surgical stent

Practice/Laboratory

When the dentist receives the stent from the laboratory, he or she can conduct a try-in to examine the aesthetics, phonetics and function of the planned dental restoration. It is advisable to make any necessary adjustments before the radiographic examination.



The radiograph serves as an additional treatment planning instrument. It enables the dental professional to establish the location and angulation of the implants, taking the bone substance and functional requirements of the future restoration into consideration.

(If necessary, radiographic balls or hollow titanium cylinders may be used with the stent.)

Application

In order to create guide channels prior to the surgical part of the treatment, the dentist (or dental technician) drills through the radiopaque teeth. These channels assist the clinician during the surgical intervention. In order to establish ideal prerequisites for the subsequent fabrication of the superstructure, the angles of the implant axes must be as parallel as possible.

CT image

Creation of guide channels

Reduced teeth as an alternative

The surgical stents are used according to the preferred method of the user with the customary materials.

Tooth and implant axes

Completed surgical stent

Implant axes on the model

Completed acrylic resin denture

Additional applications

- Insertion of metal guide channels
- Sagittal reduction of the vestibular or oral region of the guide channel up to the first half of the drill holes.

Warnings

- The inhalation of grinding dust must be avoided.
- Mucous membranes must not come into contact with uncured material.
- Prolonged or repeated skin contact with the monomer and uncured material may cause irritation and sensitivity to methacrylates.
- The Modifier Monomer contains methyl methacrylate.
- MMA is easily flammable and irritating (flash point +10 °C / 50 °F).
- Irritating to the eyes, skin and respiratory organs.
- Do not inhale vapours.
- May cause sensitization by skin contact.
- Keep away from ignition sources. Do not smoke.
- Do not empty in drains.

Storage of the Modifier Monomer / Modifier Polymer

- Store in a cool, dark place.
- Storage temperature: 12–28 °C / 54–82 °F
- See secondary packaging for storage instructions and expiry date.
- Do not use the product after the indicated date of expiration.
- Keep out of the reach of children.

Delivery form

Basic Kit

- 1x SR Vivo TAC upper anterior teeth in moulds A15, A24B, A25, A27
- 1x SR Vivo TAC lower anterior teeth in moulds A4, A8, A9
- 1x SR Ortho TAC upper posterior teeth in moulds N2, N6
- 1x SR Ortho TAC lower posterior teeth in moulds N2, N6

Modifier Kit (1 Monomer 30 ml, 1 Polymer 30 g)

Instructions for Use

Refills

- SR Vivo TAC A15, A24B, A25, A27, A4, A8, A9
- SR Ortho TAC N2, N6 upper and lower
- Polymer Modifier 30 ml
- Monomer Modifier 30 g

Source:

Pictures courtesy of M Nanni, Dr A Fanti, P Miceli, Italy

This product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent S.A.
Ivoclar Vivadent AB
Ivoclar Vivadent, Inc.