Instructions for Use

■ PRODUCT DESCRIPTION
Au-based dental metal-ceramic alloy, Type 4

■ INDICATIONS*
Inlays, Onlays, 3/4 Crowns, Crowns, Telescopic Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts / Cores, Posts, Posts with Retention, Superstructures, Partial Dentures

■ WAXING/MODELLATION
Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. The waxing design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

■ FINISHING
Provide the modeled single-tooth restoration or bridge framework with sprues of a suitable size. In general the reservoir, sprue leads, and connector sprues, whether pair shaped or traditional, must be sized according to the planned waxing design. When using the direct or indirect technique be sure that the reservoir is positioned in the heat center. The connector sprues between the reservoir and the casting should be a maximum of 2.5–3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (grams) x alloy density = grams of alloy required.

■ INVESTING
Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN-OUT
Recommended burn-out temperature: 850 °C

■ MELTING AND CASTING
Torch: Propane 0.35 bar, Oxygen 0.7 bar
Other specifics may be required by the type of casting machine. It is recommended to use a separate and separate clean burner for each metal. The recommended ratio of used material to new material is 1:1. Do not use flux.

■ CASTING TEMPERATURE
Temperature: 1260 °C

■ FRAMEWORK FINISHING
After bench cooling, clean the investment and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer for divesting. Finish the casting with carbide burs and/or with ceramic-bonded grinding instruments. Blast the surface with 50–110 micron aluminum oxide (Al₂O₃) at 2.0 bar. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.

■ OXIDATION
Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.

■ TEMPERATURE: 930 °C; **Holding time:** 5 min; **Vacuum:** No
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. In the case of Zn-containing alloys, the addition of gold to the oxide layer in 10% sulfuric acid or similar pickling solution is recommended. Note: Always observe proper safety procedures when handling acid. Before equate application clean framework carefully with water. Use the appropriate ceramic veneering material, following the manufacturer's instructions.

■ HIGHEST RECOMMENDED FIRING TEMPERATURE: 950 °C
Höchste empfohlene Brenntemperatur: 950 °C

■ HEAT TREATMENT
Hardening: 450 °C for 15 min; bench cool

■ SOLDERING AND LASER WELDING
The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool in air. Use the following soldering materials:
Pre Solder: BioPorta Lot 1020 **Flux:** High Fusing Bondal Flux
Post Solder: BioPorta OP Lot W-2 **Flux:** Bondal Flux
Laser welding Wire: BioPorta G

■ POLISHING
Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing pastes. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ CONTRAINDICATIONS
See patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician and/or allergist is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS
In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the potential side effects of dental alloys to which MRI results and to disclose the presence of dental alloys to the MRI Technician prior to conducting a test.

■ INTERACTIONS
Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION
Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and / or suitable protective masks is advised!

■ STORAGE CONDITION
Store in a dry environment at room temperature.

■ DISCLAIMER
This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products under their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

VERARBEITUNGSDATEN	
Investment material:	phosphate-bonded
Preheating/Burn-out Temperature:	850 °C
Crucible:	Carbon crucible / Ceramic crucible
Casting Temperature:	1260 °C
Oxidation:	Temperature: 930 °C; Holding time: 5 min; Vacuum: No
CTE:	(25–500 °C): 14.5 x 10 ⁻⁶ /K
Recommended Ceramic Material:	IP5 Style [®] , IP5 InLine [®]
Hardening:	450 °C for 15 min; bench cool
Pre Solder / Flux:	BioPorta Lot 1020 High Fusing Bondal Flux
Post Solder / Flux:	Porta OP Lot W-2 Bondal Flux
Laser Welding Wire:	BioPorta G

TECHNISCHE DATEN (ISO 22674:2016 & ISO 9693-1:2012)	
Type / Farbe:	4 Rich Yellow
Density (g/cm ³):	18.9
Melting Range (Solidus/Liquidus):	1040–1100 °C
Elastic Modulus (GPa):	103
Material:	Porcelain Fired
Hardened:	Hardened
Vickers Hardness:	200 210
Tickers Strength (MPa):	590 630
0.2% Proof Stress (MPa):	500 560
Elongation (%):	3 3

* See TYPE CLASSIFICATION DUE TO PHYSICAL PROPERTIES

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