Instructions for Use

EN

Gebruichsinformatioun

DE

Mode d'emploi

FR

Istruzioni d'uso

IT

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ES

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PT

Bruksanvisning

SV

Brugsanvisning

DA

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FI

■ PRODUCT DESCRIPTION

Au-based dental alloy / metal-ceramic alloy, Type 4

■ INDICATIONS*

Inlays, Onlays, Crowns, Telescopic Crowns, Conus Crowns, Bridges, Fixed Bridges, Cast Posts / Cores, Bars, Attachments, Implant Retained Superstructures, Partial Dentures

■ WAXING/MODELLATION

Design the framework in full contour or in a reduced anatomical shape taking the planned veneer into consideration. Use composite veneering material using mechanical retention. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the crown thickness design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05 – 0.2 mm.

■ SPRING

Use the included single-tooth restoration or bridge framework with springs of a suitable size. In general the reservoir, spring leads, and connector spurs, whether pair shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned so that the spring could be inserted into the reservoir and the casting should be a maximum of 2.5 – 3.0 mm in length and width. The wax pattern including the spurs must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

■ INVESTING

Use a phosphate / gypsum bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN-OUT

Recommended burn-out temperature: 700 °C

■ MELTING AND CASTING

Torch: Propan 0.15 bar, Oxygen 0.35 bar
Other specific instructions may be required by the type of casting machine. It is recommended to use a separate and clean carbon / ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended maximum safe heat for new material is 1:1. Do not use flux.
Casting Temperature: 1200 °C

■ FRAMEWORK FINISHING

After bench cooling, carefully clean and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer or other product with carbon fibers with carbon bars and/or with ceramic-bonded grinding instruments. Blast the surface with 50 – 110 micron aluminum oxide (Al₂O₃) at 2.0 bar. For veneering with Lab Composite castings, the framework must be finished and polished. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.

■ OXIDATION

Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.
Temperature: 800 °C, Holding time: 5 min; Vacuum: No
In case of oxide layer is stained, grind and blast the surface again. Repeat the oxidation in case of Zn-containing alloys. Avoid splashing of the oxide layer in 10% sulfuric acid or similar pickling solutions is recommended. Note: Always observe proper safety procedures when handling acid. Before application clean framework carefully with water. Use the appropriate ceramic cleaning material, following the manufacturer's instructions.

■ HIGHEST RECOMMENDED FIRING TEMPERATURE: 860 °C

■ HEAT TREATMENT

Hardening: 450 °C for 15 min; bench cool

■ SOLDERING AND LASER WELDING

Hardness: 15 min bei 450 °C; abkühlen lassen
The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly. Use flux sparingly.
Pre-Soldier: Porta Optimum Lot 940 Flux: High Fusing Bondal Flux
Post-Soldier: Porta Optimum Lot 710 Flux: Bondal Flux
Laser Welding Wire: Portadur IN

■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful to clean using soft brushes.

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ CONTRAINDICATIONS

For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS

In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be aware of the possibility of dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI technician prior to conducting a test.

■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and / or suitable protective masks is advised.

■ STORAGE CONDITION

Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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PROCESSING DATA	phosphate-bonded / gypsum-bonded	
Investment Material:	phosphate-bonded / gypsum-bonded	
Preheating / Burn-out Temperature:	700 °C	
Crucible:	Carbon crucible / Ceramic crucible	
Casting Temperature:	1200 °C	
Oxidation:	Temperature: 800 °C; Holding time: 5 min; Vacuum: No	
CTE:	(25 – 500 °C): 16,2 x 10 ⁻⁶ /K	
Recommended Veneering Materials:	empfohlene Verblendwerkstoffe: low-expanding, high-expanding, leucite-containing metal ceramics	
Hardening:	450 °C for 15 min; bench cool	
Pre Soldier / Flux:	Porta Optimum Lot 940 High Fusing Bondal Flux	
Post Soldier / Flux:	Porta Optimum Lot 710 Bondal Flux	
Laser Welding Wire:	Portadur IN	
TECHNISCHE DATA (ISO 22674:2016 & ISO 9693-1:2012)		
Type / Color:	4 Rich Yellow	
Density (g/cm ³):	16,0	
Melting Range (Solidus/Liquidus):	910 – 1060 °C	
Elastic Modulus (GPa):	98	
As Cast	Porcelain Fired	Hardened
Vickers Hardness:	180 190 210	210
Tickers-Härte:	180 190 210	210
Flexural Strength (MPa):	560 580 600	600
3-Point Bending Strength (MPa):	380 380 410	460
Bruchdehnung (%):	10 8 7	8
Elongation (%):	10 8 7	8

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