Description

Tetric® PowerFill is a light-curing radiopaque composite for the direct restorative treatment in posterior teeth. Tetric PowerFill cures with light in the wavelength range of 400-500 nm and can be applied in layers of up to 4 mm.

Shades

Tetric PowerFill is available in the following shades:

- Universal shades: ^{IV}A, ^{IV}B, ^{IV}W

Composition

The monomer matrix is composed of dimethacrylates (20–21 wt%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (79–80 wt.%). Additives, initiators, stabilizers and pigments are additional ingredients (< 1.0 wt%). The total content of inorganic fillers is 76–77 wt% or 53–54 vol%. The particle size of inorganic fillers is between 40 nm and 3 µm.

Indications

Conventional application (Light Intensity ≤ 2,000 mW/cm²)

- Restorations in the posterior region (Class I and II, including the replacement of individual cusps)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Reconstructive build-ups
- Restoration of deciduous teeth

Light-curing using the 3sCure mode of Bluephase® PowerCure (Light Intensity 3,050 mW/cm²)

- Restorations in the posterior region of permanent dentition (Class I and II, including the replacement of individual cusps) when light-cured from the occlusal aspect

Contraindications

The placement of Tetric PowerFill restorations is contraindicated

- If a dry working field cannot be established, or the stipulated working procedures cannot be applied.
- If a patient is known to be allergic to any of the ingredients of Tetric PowerFill

Safety notes

- Do not place light in direct contact with unprotected gingiva, mucous membrane or skin.
- The 3sCure curing mode must not be used in case of caries profunda and very deep cavities.

Side effects

In individual cases, components of Tetric PowerFill may lead to sensitization. Tetric PowerFill should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp and cover with suitable cavity liner).

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Interactions

Materials containing phenolic substances, e.g. eugenol/clove oil, inhibit the polymerization of methacrylate-based materials. Consequently, application of such materials together with Tetric PowerFill must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application

1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

Cavity preparation is carried out according to the requirements of the adhesive technique, i.e. protecting the tooth structure. Do not prepare sharp, internal edges and angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. For the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

The 3sCure mode must not be used for restorations in areas close to the pulp. Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use a pressure-resistant cement (e.g. a glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® Universal, a single-component, light-cured adhesives for direct and indirect bonding procedures that features compatibility with all etching techniques.

7. Application of Tetric PowerFill

- In order to achieve optimum results, Tetric PowerFill should be applied in increments of max. 4 mm and adapted to the cavity walls with suitable instrument (e.g. OptraSculpt®).
- Sufficient exposure to the curing light prevents incomplete polymerization.
- For the recommendations regarding exposure time per increment and light intensity see table 1.

The instructions for use of the curing light must be observed.

- When using a metal matrix, additionally polymerize the composite material from the buccal or lingual/palatal aspect after removing the matrix, if no Bluephase® curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a diverging scattering angle, the composite material should also be light-cured again (no 3sCure mode), but observing the above mentioned safety notes.
- As an option, a flowable composite (e.g. Tetric EvoFlow®, Tetric EvoFlow® Bulk Fill or Tetric PowerFlow) can be applied. This optional layer must be cured separately (please refer to the respective Instructions for Use).

8. Finishing / Checking the occlusion / Polishing

Remove excess material with suitable finishers or fine diamonds after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric PowerFill can be used in combination with Tetric PowerFlow, Tetric EvoCeram, Tetric EvoFlow and Tetric EvoFlow Bulk Fill. Use the 3sCure mode of Bluephase PowerCure for Tetric PowerFill or Tetric PowerFlow materials only.
2. In the case of repairs, additional Tetric PowerFill can be directly applied to the polymerized material. If the Tetric PowerFill restoration has already been polished, it must first be roughened using diamond grinding instrument and wetted with Adhese Universal before fresh Tetric PowerFill is applied.
3. Tetric PowerFill should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If Tetric PowerFill is applied from the Cavifil directly into the oral cavity of the patient, the Cavifil should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
5. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.
6. The recommended increment thickness is based on hardness profile measurements.

Warning

- Unpolymerized Tetric PowerFill should not come in contact with skin, mucous membrane and eyes. Unpolymerized Tetric PowerFill may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Storage temperature 2 – 28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use Tetric PowerFill after the indicated date of expiration.
- Date of expiration: see information on Cavifils, syringes and packages.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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