

The Heliomolar[®] Family:

Heliomolar[®]

Heliomolar[®] Flow

Heliomolar[®] HB

Scientific Documentation

Table of Contents

1. Introduction	3
1.1 Heliomolar.....	3
1.2 Heliomolar Flow	3
1.3 Heliomolar HB.....	3
2. Technical Data	4
3. Physical Studies	5
3.1 Abrasion.....	5
3.2 Shade stability.....	5
3.3 Physical data	6
3.4 Fluoride release.....	6
3.5 Polymerization shrinkage.....	6
3.6 Marginal adaptation	7
4. Clinical Studies	8
4.1 Heliomolar.....	8
4.2 Heliomolar HB.....	9
5. Toxicology	11
5.1 Acute oral toxicity	11
5.2 Compatibility with the mucous membrane and local irritation of skin	11
5.3 Elution tests.....	11
5.4 Ytterbium trifluoride.....	11
5.5 Sensitization	11
5.6 Mutagenic properties.....	11
5.7 Literature on toxicology	12
6. Literature	13

1. Introduction

In the last few years, the growing call for invisible restorations and the search for alternative materials to amalgam have led to an increase in the demand for composite materials. In 1984, Ivoclar Vivadent introduced Heliomolar, which has become one of the most widely accepted dental composite materials. In 1994, CRA wrote that Heliomolar was being successfully used all over the world since its introduction in 1984 and even eleven years after its launch it continued to be one of the best composite materials. In April 1999, Heliomolar Flow, a flowable version of Heliomolar, was introduced to cover all the requirements of dentists in terms of indications and application procedures. Furthermore, in December 2000, a condensible product variant called Heliomolar HB was launched, completing the range of Heliomolar materials.

1.1 *Heliomolar*

Heliomolar falls into the category of inhomogeneous microfilled composites. In dentistry, microfillers are materials whose filler particles are smaller than 1µm. Prepolymers may be added to the microfilled composite to enhance its consistency and physical properties as well as to increase its filler content. The prepolymers used in Heliomolar are microfilled pre-polymerized composites that exhibit virtually the same properties as the matrix.

1.2 *Heliomolar Flow*

Heliomolar Flow is a flowable version of Heliomolar. The monomer content is slightly higher than that of the original Heliomolar to render the material flowable. Because of its flowable consistency, Heliomolar Flow is particularly indicated for Class V defects, mini-cavities of all kinds, preventive resin restorations as well as the repair of composite and ceramic veneers. Furthermore, many dentists use Heliomolar flow as a first thin increment under Heliomolar or Heliomolar HB restorations because the flowable consistency facilitates adaptation to the cavity bottom and walls.

1.3 *Heliomolar HB*

Heliomolar HB is the latest member of the Heliomolar family. The adjunct HB stands for Heavy Body. Heliomolar HB falls into the category of what are known as packable or condensible composites. Ivoclar Vivadent offers Heliomolar in three different consistencies to meet the varying requirements placed on the handling properties of composites. Heliomolar HB is particularly suitable for direct restorations in the posterior region.

Two slight alterations of the original Heliomolar were required to obtain the heavy-body consistency of Heliomolar HB. Firstly, the material was rendered less sticky by slightly modifying the proportional composition of the monomer mixture, ie the portion of the comparatively large copolymers was lowered while the portion of microfillers was slightly raised. As a result, the viscosity of the material increased. Secondly, a rheology modifier in the form of an organically modified compound silicate was added to the material. The compound silicate, which contains surface linked, long chain organic groups, increases the firmness of Heliomolar HB but does not compromise the material's modelling properties.

2. Technical Data

Standard composition	Heliomolar	Heliomolar HB	Heliomolar Flow
Bis-GMA, Urethane dimethacrylate	19	18	32
Decandiol dimethacrylate	3	5	-
Triethylene glycol dimethacrylate	-	-	8
Highly dispersed silicon dioxide, Prepolymer, Ytterbium trifluoride	77	76	59
Stabilizers, catalysts and pigments	< 1	< 1	< 1

(Figures in wt%)


Physical properties	Unit	Heliomolar	Heliomolar HB	Heliomolar Flow
Flexural strength	MPa	100	125	96
Modulus of elasticity	MPa	6000	6500	4400
Compressive strength	MPa	340	315	260
Vickers hardness	MPa	350	415	300
Water absorption	$\mu\text{g}/\text{mm}^3$	25	24.9	< 30
Water solubility	$\mu\text{g}/\text{mm}^3$	1	1.5	< 3
Radiopacity	% Al	250	265	200
Total filler content	wt %	66.7	66.7	51
Total filler content	vol %	46	46	30

In accordance with: ISO 4049 – Polymer-based filling, restorative and luting materials


3. Physical Studies

3.1 Abrasion


When Heliomolar was introduced, it set a new standard in terms of resistance to abrasion of composite materials. In 1987, Tani et al published a study, in which Heliomolar’s high resistance to abrasion was confirmed. Furthermore, in a study conducted by Leinfelder in 1991, Heliomolar was shown to be one of the materials that exhibited the lowest wear. The experimental of this study included 400,000 cycles of repeated stress and PMMA pellets as the abrasive. While the majority of the composite materials against which Heliomolar was compared in 1987 and 1991 are no longer available, Heliomolar as well as Heliomolar Flow and Heliomolar HB, the two new Heliomolar versions, still rank among the leading composites in terms of resistance to abrasion (Sorensen 2000, Report on file).



Tani et al (1987)



Leinfelder (1991)



Sorensen (2000)

Wear of Heliomolar, measured in three independent studies: The results of Tani and Leinfelder are taken from the above-mentioned studies. The results of Sorensen are on file.

3.2 Shade stability

Heliomolar exhibited a high degree of shade stability in trials involving artificial ageing (Powers et al, 1988). With regard to most of the parameters examined, Heliomolar demonstrated the least change when subjected to artificial ageing.

	ΔY	ΔS	ΔK	ΔCR
Heliomolar	-1.6	-0.02	0.021	0.003
Herculite	-1.9	0.13	0.043	0.070
P-30	-9.3	0.07	0.017	0.068
Bis-Fil II	-2.3	0.04	0.039	0.021
Estilux Posterior	4.8	0.32	-0.006	0.122
Distalite	1.6	0.03	-0.14	0.006

Optical stability of Heliomolar

ΔY = Change in the degree of transmission
 ΔK = Change in the absorption coefficient

ΔS = Change in light scattering
 ΔCR = Change in opacity


3.3 Physical data

	Compressive strength [N/mm ²]	Flexural strength [N/mm ²]	Flexural modulus [GPa]
Heliomolar HB	380	114	7.3
Tetric Ceram HB	338	138	13.2
P-60	524	170	15
Solitaire	412	54	3
Alert	399	119	10
SureFil	426	131	11
Prodigy Condense	371	104	8
Tooth enamel	384	90	84

Munoz (Loma Linda University, California, USA)

3.4 Fluoride release

Arends und Ruben (1988) measured the fluoride releasing capabilities of dental composites. Heliomolar was proven to continuously release fluoride for over one year. Moreover, Arends and Zee (1990), who conducted a study using an artificial mouth model, showed that dentin and enamel take up the fluoride released by Heliomolar



Fluoride release of Heliomolar (Arends and Ruben 1988)

3.5 Polymerization shrinkage

Among all the composites commercially available, Heliomolar is one of the materials that demonstrates the lowest polymerization shrinkage (Feilzer et al, 1988; Soltez, 2000).


	Feilzer et al, 1988	Soltez, 2000
Heliomolar	2.2 ± 0.1	1.96 ± 0.07
Heliomolar HB		2.19 ± 0.10
Tetric Ceram		2.76 ± 0.05
Tetric Ceram HB		2.55 ± 0.01
Esthet X		2.62 ± 0.04
Point 4		3.01 ± 0.10

Brilliant Lux	3.5 ± 0.0	
Clearfil Posterior	4.5 ± 0.5	
Herculite	3.0 ± 0.2	
P-30	2.6 ± 0.3	
Prisma Fil	3.3 ± 0.2	

Polymerization shrinkage in vol %, 1 hour after curing. The data of Soltez (Fraunhofer Insitut Werkstoffmechanik, Freiburg, Germany) were measured for Ivoclar Vivadent.

3.6 Marginal adaptation

In a study conducted by Stähle and Ackermann in 1991, 70 percent of the occlusal restorations fabricated of Heliomolar demonstrated tight margins after they were exposed to cyclic loading at 200 N. By contrast, none of the gold inlays, which were inserted using zinc phosphate, exhibited tight margins.



Depth of dye penetration:

0 = no dye penetration

1 = dye penetration into the upper half of the restoration

2 = dye penetration into the lower half of the restoration

3 = dye penetration down to the cavity floor

Study on marginal integrity by Stähle and Ackermann (1991)

4. Clinical Studies

4.1 Heliomolar

Heliomolar was launched in 1984. Since then, several three- and five-year studies on the material have been accomplished.

3-year studies

Knibbs and Smart (1992) conducted a comparative three-year study with Heliomolar and amalgam. Fifty-two fillings of each material were placed and monitored in terms of marginal adaptation, surface quality, anatomical shape, and proximal antagonist contact. No significant difference in the performance of Heliomolar and amalgam was found; both materials produced favourable clinical results.

Lundin et al (1990) monitored the clinical performance of six different dental composites over three years. The materials, which were placed in Class II cavities, showed very low failure rates. Heliomolar demonstrated the lowest abrasion of all six materials

Taylor et al (1994) summarized the data of ten clinical three-year studies on the wear of restorative materials. Heliomolar ranked second among the 28 composites examined.

In a clinical three-year study conducted by CRA in the US, Heliomolar produced the best results among the 21 restorations investigated. The following aspects were monitored: abrasion, marginal adaptation, surface quality, antagonist abrasion, resistance to fracture and shade adaptation.

Clinical Research Associates Newsletter, Volume 18, Issue 5 May 1994: Comparative performance of 21 class 2 materials at 3 years

5-year studies

In a clinical study conducted by Leinfelder (Mazer and Leinfelder, 1992), 68 Class I and II restorations were placed using Heliomolar. The mean abrasion was as low as 7.7 µm/year.. According to Leinfelder, the Heliomolar restorations were characterized by a high degree of shade stability, surface smoothness, resistance to abrasion and acceptance by patients throughout the five years. During the first two years, none of the restorations demonstrated secondary caries. Only two restorations showed signs of incipient secondary caries. Leinfelder attributes the low rate of secondary caries to the fact that Heliomolar releases fluoride ions.

	Abrasion after 2 years	Secondary caries after 2 years
Heliomolar	12 µm	0 %
Herculite	30 µm	3 %
Bisfil-1	42 µm	4.5 %
P-10	135 µm	2 %

The authors compared the data on abrasion and secondary caries with the corresponding data of other materials. Heliomolar restorations demonstrated the lowest abrasion and were the least prone to developing secondary caries (Mazer and Leinfelder, 1992).

Setcos and Phillips (1995) compared Heliomolar and amalgam in a clinical study on Class I and II restorations. The observation period was five years and the following aspects were monitored: shade stability, marginal discoloration, anatomic shape, marginal integrity, secondary caries, abrasion and sensitivity to temperature. The results prove the favourable clinical performance of Heliomolar.

A CRA Newsletters (1989) summarized the results of clinical trials conducted on 21 different dental composites placed in Class II cavities. In the trials, which lasted between two to five years, the following aspects were monitored: marginal adaptation, proximal contact areas, postoperative sensitivity, secondary caries, surface quality, discoloration and shade. The Heliomolar restorations were attested favourable ratings with regard to all of these aspects. Moreover, Heliomolar demonstrated the lowest abrasion of all the materials tested.

If the results of the 3- and 5-year studies are all added up and a loss analysis according to Kaplan-Meier is conducted, very high survival rates are observed.

	3-year studies	5-year studies
Base line	100%	100%
1 year	99.6 ± 0.4	100%
2 years	98.0 ± 0.9	100%
3 years	96.3 ± 1.2	98 ± 1%
4 years		97 ± 2%
5 years		93 ± 2%

Conclusion: Heliomolar is one of the most proven and clinically successful dental composites on the market.

4.2 Heliomolar HB

Clinical studies to examine Heliomolar HB have been initiated. The 12-month results are available.

Head of study: Dr Jim R Dunn, D Carlos Munoz
Loma Linda University, California, USA

Objective: Examine the clinical performance of Heliomolar HB

Experimental: Fifty Class II cavities were placed in a total of 32 patients, using Heliomolar HB and Excite dentin adhesive. Heliomolar HB was applied in increments. The restorations were examined after six and twelve months. The next follow-up evaluation is due to be conducted 24 months after placement.

Results:

Criteria	Heliomolar HB		
	1 months	6 months	12 months
Restorations evaluated	44	48	43
Marginal adaptation	100% A	100% A	88% A
Interproximal anatomical shape	98% A	100% A	95% A
Post-operative sensitivity	90% A	96% A	91% A
Secondary caries	100% A	100% A	100% A
Marginal discoloration	96% A	98% A	88% A
Proximal contact areas	86% A	96% A	91% A
Surface polish	100% A	100% A	100% A
Retention	100 %	100%	98% A

Conclusion: The evaluation conducted six months after placement showed excellent results for Heliomolar HB. After 12 months only one restoration had to be replaced and no other Charlie ratings were observed.

Head of study: Dr Giovanni Dondi dall'Orologio
University of Bologna, Italy

Objective: Examine the clinical performance of Heliomolar HB, placed in Class I and II cavities. Place an initial base layer of Heliomolar Flow in half of all the restorations. After 6 months all patients and after 12 months 62 patients could be evaluated. The next recall is planned after 24 months.

Experimental: Restorations fabricated of Heliomolar HB and Heliomolar Flow/Heliomolar HB were placed in 62 patients according to a split mouth design. Excite was used as the adhesive. Heliomolar HB was placed in increments. The restorations were evaluated after six months. The next follow-up evaluation is due to be conducted at 12 months after placement.

Results:

Criteria	Heliomolar HB		
	Baseline	6 months	12 months
Marginal integrity	100% A	100% A	87% A, 13% B
Discoloration	100% A	100% A	100% A
Secondary caries	100% A	100% A	100% A
Surface quality	100% A	100% A	100% A
Anatomical shape	100% A	100% A	100% A
Post-operative sensitivity	97% A, 3% B	97% A, 3% B	91% A, 9% B
Retention	100% A	100% A	100% A

Criteria	Heliomolar Flow & Heliomolar HB		
	Baseline	6 months	12 months
Marginal integrity	100% A	100% A	100% A
Discoloration	100% A	100% A	100% A
Secondary caries	100% A	100% A	100% A
Surface quality	100% A	100% A	100% A
Anatomical shape	100% A	100% A	100% A
Post-operative sensitivity	98% A, 2% B	98% A, 2% B	100% A
Retention	100% A	100% A	100% A

Conclusion: The evaluation conducted six months after placement shows favourable results for both kinds of Heliomolar HB restorations, i.e., those placed with Heliomolar Flow and those placed without Heliomolar Flow. The 12 months results confirm the excellent clinical success of Heliomolar HB. The data indicate that in the long term the use of Heliomolar Flow as first thin increment may result in better margins.

5. Toxicology

Heliomolar is applied directly into the cavity. The resin is light-cured in the oral cavity.

The following toxicological risks had to be investigated:

- Acute oral risk: patients may accidentally swallow the portion of uncured Heliomolar placed in the cavity
- Local incompatibility with surrounding tissue that comes into contact with the material
- Possible sensitizing reactions
- Peroral long-term risk by eluted low-molecular components
- Mutagenic potential of eluted low-molecular components

5.1 Acute oral toxicity

Tests on the acute oral toxicity of the uncured formulation produced the following values:

LD₅₀ (rats) oral p.o. > 5000 mg/kg [1]

ISO filler RO = HELIOMOLAR RO without ytterbium trifluoride

LD₅₀ (rats) p.o. > 5000 mg/kg [2]

Ytterbium trifluoride

LD₅₀ (rats) p.o. > 5000 mg/kg [3]

5.2 Compatibility with the mucous membrane and local irritation of skin

A primary irritation index of 0.3 was measured in the mucuous membrane. On the basis of this result, the material may be regarded as having a minimally irritating effect in an uncured state [4]. The skin irritation index was also 0.3, indicating a minimal primary irritation effect on skin [5].

5.3 Elution tests

A mutation assay was conducted to assess the effect of oral long-term exposure to possible eluates [6, 7]. Total migration of 230 µg/cm² was measured.

5.4 Ytterbium trifluoride

The toxicity of ytterbium trifluoride in Heliomolar RO is discussed in a comprehensive report by Dr Manfred Herbst [8]. This report reaches the conclusion that ytterbium trifluoride does not involve any health risk if it is used as a component of a polymethacrylate-based filling material.

5.5 Sensitization

Uncured Heliomolar was subjected to a sensitization (maximization) test in guinea pigs. In this rigorous test [9], slight irritation occurred. However, it may be assumed that the cured material has no sensitizing effect.

5.6 Mutagenic properties

Mutation did not occur in an Ames test (reversible mutagenicity assay) under the conditions chosen, neither in the HGPRT genes of V79 cells (Chinese hamster) [10] nor in Salmonella typhimurium strains (TA 1535, TA 1537, TA 1538, TA 98 and TA 100) [11]. Furthermore, in a

chromosome aberration assay conducted on Chinese hamsters, no chromosomal mutations occurred in correlation with Heliomolar [12]. On the basis of these tests, mutagenic effects of Heliomolar are not indicated.

Conclusion:

On the basis of the data on hand, an acute or chronic risk for patients is not indicated, if Heliomolar is used properly.

The results of the clinical studies and the general toxicological data on composite materials support the above conclusion.

The chemistry of Heliomolar Flow and Heliomolar HB are similar to that of Heliomolar. The monomer content and the filler composition have only been slightly changed to modify the consistency of the materials.

Thus, the toxicological data on Heliomolar also apply to Heliomolar Flow and Heliomolar HB.

5.7 Literature on toxicology

- | | |
|---|---|
| [1] RCC Project 048857: Acute Oral Toxicity (LD ₅₀) Study with Heliomolar Radiopaque in rats. July 1, 1985. | Zahnfüllungen, Swiss Dent 12, 15-19, 1984 |
| [2] RCC Project 034593: Acute Oral Toxicity (LD ₅₀) Study with Isofüller RO in rats. August 31, 1984. | [8] Dr. Manfred Herbst: Assessment: Toxicity of Ytterbiumtrifluoride in Heliomolar Radiopaque. Liestal, July 20, 1985 |
| [3] RCC Project 048881: Acute Oral Toxicity (LD ₅₀) Study with ytterbium trifluoride, anhydrous in rats. July 1, 1985 | [9] RCC Project 205323: Contact Hypersensitivity to Heliomolar Radiopaque in Albino Guinea Pigs: Maximisation Test. April 21 1988 |
| [4] RCC Project 048868: Primary Eye Irritation Study with Heliomolar Radiopaque in rabbits. July 30, 1985 | [10] CCR Project 296324: Salmonella Thyphimurium Reverse Mutation Assay with Heliomolar Radiopaque Universal. August 1992 |
| [5] RCC Project 048870: Primary Skin Irritation Study with Heliomolar Radiopaque in rabbits. July 30, 1985 | [11] CCR Project 313918: Gene Mutation Assay in Chinese Hamster V79 Cells in vitro with Heliomolar Radiopaque Universal. February 18, 1993 |
| [6] Prüfbericht: Migration von Heliomolar und HELIOMOLAR RO nach AP-22. A. Weber, 11. März 1985, Schaan FL | [12] CCR Project 313920: Chromosome Aberration Assay in Chinese Hamster V79 Cells in vitro with Heliomolar Radiopaque Universal. April 26, 1993 |
| [7] Leimgruber R, Thöny D, Weber A: Elutionstest zur Prüfung auf Anwesenheit extrahier-barer Bestandteile von | |

6. Literature

Arend J, van der Zee Y

Fluoride uptake in bovine enamel and dentin from fluoride-releasing composite resin
 Quintessence International 21 (1990) 541-544

Arends J, Ruben J

Fluoride release from a composite
 Quintessence International 19 (1988) 513-514

Clinical Research Associates

Comparative performance of 21 class 2 materials at 3 years
 Clinical Research Associates Newsletter, Volume 18, Issue 5 May 1994

Clinical Research Associates

Restorative resins, new products, one year clinical performance
 Clinical Research Associates Newsletter, Volume 18, Issue 11 November 1994

Feilzer AJ, De Gee AJ, Davidson CL

Curing contraction of composites and glass-ionomer cements
 J Prosthetic Dent 59 (1988) 297-300

Knibbs PJ, Smart ER

The clinical performance of a posterior composite resin restorative material, HELIOMOLAR RO: 3-year report
 J Oral Rehabil 19(1992) 231--237

Leinfelder KF, Mirshaahidi M, Cury C, O'Neal W

An in vitro wear device for determining wear of posterior composites
 J Dent Res 70 (1991) 345

Publications on Heliomolar and Heliomolar Flow

Restorative resins, new products, status report #1 -clinical characteristics
 CRA Newsletter 17 (1993) 0

Alternatives for class 3 restorations (results of CRA clinical trial of 21 materials at 3 years)
 CRA Newsletter 18 (1994) 1-4

Alhadainy HA, Abdalla AI

2-year clinical evaluation of dentin bonding systems
 Am J Dent 9 (1996) 77-79

Arends J, Ruben J, Dijkman AG

The effect of fluoride release from a fluoride-containing composite resin on secondary caries: an in vitro study
 Quintessence Int 21 (1990) 671-674

Arends J, Van der Zee Y

Fluoride uptake in bovine enamel and dentin from a fluoride-releasing composite resin
 Quintessence Int 21 (1990) 541-544

Ludin SA, Andersson B, Koch G, Rasmusson CG
Class II composite resin restorations: A three year clinical study of six different posterior composites

Swed Dent J 14 (1990) 104-115

Mazer R, Leinfelder K

Evaluating a microfill posterior composite resin - a five year study
 J Am Dent Assoc 123 (1992) 33-38

Powers JM, Bakus ER, Goldberg AJ

In vitro color changes of posterior composites
 Dent Mater 4 (1988) 151-154

Setcos JC

Heliomolar radiopaque als Amalgamersatz?
 Phillip J 12 (1995) 93-95

Staehele HJ, Ackermann J

Experimentelle Untersuchungen über die Dichtigkeit von okklusalen Kompositfüllungen im Vergleich zu okklusalen Metallinlays
 ZWR 100 (1991) 648-653

Tani Y, Goto H, Ida K

Wear of posterior composite resins
 Dental materials Journal 6 (1987) 165-174

Taylor DF, Bayne SC, Leinfelder KF, Davis S, Koch GG

Pooling of long term clinical data for posterior composites
 Am J Dent 7 (1994) 167-174

Ashe MJ, Tripp GA, Eichmiller FC, George LA, Meiers JC

Surface roughness of glass-ceramic insert - composite restorations: assessing several polishing techniques
 J Am Dent Assoc 127 (1996) 1495-1500

Bayne SC, Heymann HO, Edward J, Swift JR

Update on dental composite restorations
 J Am Dent Assoc 125 (1994) 687-701

Bollen CML, Lambrechts P, Quiryen M

Comparison of surface roughness of oral hard materials to the threshold surface roughness for bacterial plaque retention: A review of the literature
 Dent Mater 13 (1997) 258-269

Bose M, Ott KHR

Glättung von (Füllungs-) Werkstoffen, Zahnschmelz und Dentin durch Prophylaxepasten in vitro
 Dtsch Zahnärztl Z 50 (1995) 840-843

- Bouschlicher MR, Cobb DS, Boyer DB
Radiopacity of compomers, flowable and conventional resin composites for posterior restorations
Oper Dent 24 (1999) 20-25
- Bryant RW, Hodge KLV
A clinical evaluation of posterior composite resin restorations
Aust Dent J 39 (1994) 77-81
- Bryant RW, Marzbani N, Hodge KLV
Occlusal margin defects around different types of composite resin restorations in posterior teeth
Oper Dent 17 (1992) 215-221
- Castline JM, TU SJ
Measured light intensities at simulated dentin/composite bonding interfaces
J Dent Res 79 (2000) 280
- Christensen GJ
Restorative Dentistry: An update for Practitioners, Educators, Examining Boards
J Am Dent Assoc 126 (1995) 1165-1168
- Christensen GJ, Christensen RP
A new technique for the restoration of worn anterior teeth - 1995
J Am Dent Assoc 126 (1995) 1543-1546
- Chung K
Effects of finishing and polishing procedures on surface texture of resin composites
Dent Mater 10 (1994) 325-330
- Condon JR, Ferracane JL
Assessing the effect of composite formulation on polymerization stress
J Am Dent Assoc 131 (2000) 497-503
- Condon JR, Ferracane JL
Polymerization contraction stress of commercial composites
J Dent Res 77 (1998) 639
- Condon JR, Ferracane JL
Evaluation of composite wear with a new multi-mode oral wear simulator
Dent Mater 12 (1996) 218-226
- Davidson CL, Abdalla AI
Effect of thermal and mechanical load cycling on the marginal integrity of Class II resin composite restorations
Am J Dent 6 (1993) 39-42
- De Backer J, Dermaut L
Visible light sources and posterior visible light cured resins: a practical mixture
Quintessence Int 17 (1986) 635-641
- De Goes MF, Garcia-Godoy F, Cardenas L
Linear and volumetric polymerization shrinkage of the restorative materials
J Dent Res 78 (1999) 484
- Dietschi D, Ciucchi B, Holz J
A clinical trial of four light curing posterior composite resins: 9-month report
Quintessence Int 20 (1989) 641-652
- Dietschi D, Holz J
Klinische Untersuchung von vier lichthärtenden Kompositmaterialien für den Seitenzahnbereich - Ergebnisse nach zwei Jahren (1)
Quintessenz 42 (1991) 743-751
- Dietschi D, Magne P, Holz J
Recent trends in esthetic restorations for posterior teeth
Quintessence Int 25 (1994) 659-677
- Dijkman GEHM, Arends J
Secondary Caries in situ around Fluoride-Releasing Light-Curing Composites: A Quantitative Model Investigation on Four Materials
Caries Res 26 (1992) 351-357
- Dijkman GEHM, De Vries J, Lodding A, Arends J
Long-term fluoride release of visible light-activated composites in vitro: a correlation with in situ demineralisation data
Caries Res 27 (1993) 117-123
- Donly K, Gomez C
In vitro demineralization-remineralization of enamel caries at restoration margins utilizing fluoride-releasing composite resin
Quintessence Int 25 (1994) 355-358
- Fallo GJ, Wakefield CW
Effects of uncontrolled outdoor storage on the polymerization, manipulation, and appearance of visible light-cured composite resin and resin-modified glass ionomer materials
Military Medicine 161 (1996) 290-293
- Feilzer AJ, De Gee AJ, Davidson CL
Relaxation of Polymerization Contraction Shear Stress by Hygroscopic Expansion
J Dent Res 69 (1990) 36-39
- Feilzer AJ, De Gee AJ, Davidson CL
Curing contraction of composites and glass-ionomer cements
J Prosthet Dent 59 (1988) 297-300
- Ferracane JL
Current Trends in Dental Composites
Crit Rev Oral Biol Med 6 (1995) 302-318
- Ferracane JL, Condon JR
Post-cure heat treatments for composites: properties and fractography
Dent Mater 8 (1992) 290-295
- Ferracane JL, Mitchem JC, Condon JR, Todd R
Wear and marginal breakdown of composites with various degrees of cure
J Dent Res 76 (1997) 1508-1516
- Francci C, Perdigo J, Cardoso PEC, Meira JBC, Nunes MF

The effect of composite resin, adhesive system, and low-viscosity liner on microleakage

J Dent Res 78 (1999) 388

Freilich MA, Goldberg AJ, Gilpatrick RO, Simonsen RJ

Three-year occlusal wear of posterior composite restorations

Dent Mater 8 (1992) 224-228

Frenzel C, Viohl J

Biegefestigkeit von neun Kompositen in Abhängigkeit von verschiedenen Reparaturverfahren

Dtsch Zahnärztl Z 49 (1994) 729-732

Fujimitsu T, Kato Htoh K, Wakumoto S

The adaptation of composite to the dentin cavity wall

Dent Mater J 8 (1989) 141-146

Fukushima M, Setcos JC, Phillips RW

Marginal fracture of posterior composite resins

J Am Dent Assoc 117 (1988) 577-583

Fuqua S, Burgess JO, Chan DCN, Nunez A

Recurrent caries in OF restorations in a rodent model

J Dent Res 77 (1998) 118

Gogswaardt DC

Die Polymerisationstemperatur lichtvernetzbarer Seitenzahn-Kompositmaterialien in Abhängigkeit von der Schichtstärke

Zahnärztl Praxis 5 (1995) 10-12

Haak R, Noack MJ

Möglichkeiten der Röntgendiagnostik bei Amalgamersatzmaterialien

Quintessenz 47 (1996) 1551-1559

Hellwig E, Klimek J, Achenbach K

Auswirkung der Schichttechnik auf die Polymerisation von zwei lichthärtenden Komposit-Füllungsmaterialien

Dtsch Zahnärztl Z 46 (1991) 270-273

Herrgott AM, Ziemiecki TL, Dennison JB

An evaluation of different composite resin systems finished with various abrasives

J Am Dent Assoc 119 (1989) 729-732

Hibst R, Keller U

Randschluss und Haftfestigkeit von Kompositfüllungen - Präparation und Konditionierung mit dem Er:YAG-Laser

ZWR 104 (1995) 78-83

Hickel R, Kunzelmann KH, Obermeier T

Die Kompositfüllung im Seitenzahnbereich - Teil I & Teil II

ZWR 103 (1994) 610-618

Kanca III J

Maximizing the cure of posterior light-activated resins

Quintessence Int 17 (1986) 25-27

Kanca III J

Visible light-activated composite resins for posterior use - a comparison of surface hardness and uniformity of cure. Update

Quintessence Int 16 (1985) 687-690

Kawai K, Leinfelder K

In vitro evaluation of OCA wear resistance of posterior composites

Dent Mater 11 (1995) 246-251

Knibbs PJ, Smart ER

The clinical performance of a posterior composite resin restorative material, Heliomolar RO: 3-year report

J Oral Rehab 19 (1992) 231-237

Krejci I, Stergiou G, Lutz F

Einfluss der Nachvergütung auf die Verschleissfestigkeit von Kompositmaterialien

Dtsch Zahnärztl Z 46 (1991) 400-406

Lambert D

Closing diastemas, replacing small amalgams using Heliomolar einforced microfill restoratives

Dental Products Report 0 (2000) 92-93

Leinfelder K

Posterior composite resins: The materials and their clinical performance

J Am Dent Assoc 126 (1995) 663-676

Leinfelder K

Using composite resin as a posterior restorative material

J Am Dent Assoc 122 (1991) 65-70

Leinfelder K, Beaudreau RW, Mazer R

An in vitro device for predicting clinical wear

Quintessence Int 20 (1989) 755-761

Liebenberg WH

Direct access to equivocal approximal carious lesions

Quintessence Int 27 (1996) 607-617

Lundin SA, Andersson B, Koch G, Rasmusson CG

Class II composite resin restorations: a three-year clinical study of six different posterior composites

Swed Dent J 14 (1990) 105-114

Mair LH

The silver sorption layer in dental composites: three year results

Dent Mater J 15 (1999) 408-412

Mair LH

An investigation into the permeability of composite materials using silver nitrate

Dent Mater 5 (1987) 109-114

Mair LH

Subsurface compression fatigue in seven dental composites

Dent Mater 10 (1994) 111-115

Mair LH, Vowles RW

The effect of thermal cycling on the fracture toughness of seven composite restorative materials

Dent Mater 5 (1989) 23-26

Margolis HC, Moreno EC, Murphy BJ

Effect of low levels of fluoride in solution on enamel demineralization in vitro

J Dent Res 65 (1986) 23-29

Mazer R, Leinfelder K

Evaluating a microfill posterior composite resin - a five year study

J Am Dent Assoc 123 (1992) 33-38

Mazer R, Leinfelder K, Russell CM

Degradation of microfilled posterior composite

Dent Mater 8 (1992) 185-189

Milosevic A

The influence of surface finish and in-vitro pellicle on contact-angle measurement and surface morphology of three commercially available composite restoratives

J Oral Rehab 19 (1992) 85-97

Momoi Y, McCabe JF

Hygroscopic expansion of resin based composites during 6 months of water storage

Br Dent J 176 (1994) 91-96

Mount GJ

The tensile strength of the union between various glass ionomer cements and various composite resins

Aust Dent J 34 (1989) 136-146

Munksgaard EC, Nolte J, Kristensen K

Adherence of chewing gum to dental restorative materials

Am J Dent 8 (1995) 137-139

NaBadalung DP, Nicholls JI, Brudvik JS

Frictional resistance of removable partial dentures with retrofitted resin composite guide planes

Int J Prosthodont 10 (1997) 116-122

Oden A, Ruyter IE, Oysaed H

Creep and recovery of composites for use in posterior teeth during static and dynamic compression

Dent Mater 4 (1988) 147-150

Palenik CJ, Setcos JC

Antimicrobial abilities of various dentine bonding agents and restorative materials

J Dent 24 (1996) 289-295

Papagiannoulis L, Tzoutzas J, Eliades G

Effect of topical fluoride agents on the morphologic characteristics and composition of resin composite restorative materials

J Prosthodont 77 (1997) 405-413

Pelka M, Ebert J, Schneider H, Krämer N, Petschelt A

Comparison of two- and three-body wear of glass-ionomers and composites

Eur J Oral Sci 104 (1996) 132-137

Pelka M, Frankenberger R, Sindlinger R, Petschelt A

Verschleißsimulation natürlicher Zahnhartsubstanzen im abrasiven Kontakt

Dtsch Zahnärztl Z 53 (1998) 61-64

Pilo R, Brosh T, Lugassy H, Baharav H, Helft M

The effect of irradiation time on the shear strength of composites

Dent Mater 10 (1994) 338-342

Powers JM, Bakus ER, Goldberg AJ

In vitro color changes of posterior composites

Dent Mater 4 (1988) 151-154

Powers JM, McDowell GC, Lang BR

In vivo wear. Part II: Wear and abrasion of composite restorative materials

J Prosthodont 60 (1988) 242-249

Prati C

Early marginal microleakage in class II resin composite restorations

Dent Mater 5 (1989) 392-398

Rasmusson CG, Lundin SA

Class II restorations in six different posterior composite resins: Five-year results.

Swed Dent J 19 (1995) 173-182

Rees JS, Jacobsen PH

The polymerization shrinkage of composite resins

Dent Mater 5 (1989) 41-44

Reich E, Schmalz G, Federlin M

Randspaltverhalten von Keramik- und Kompositinlays in vitro

Dtsch Zahnärztl Z 45 (1990) 656-660

Reinhardt KJ

Der Einfluss der Lichtquelle auf die Randständigkeit von Kompositfüllungen

Dtsch Zahnärztl Z 46 (1991) 132-134

Reinhardt KJ

Der Einfluss des Bonders auf die Randständigkeit und Belastbarkeit von Komposit-Füllungen

Dtsch Zahnärztl Z 47 (1992) 176-178

Retief DH, Mandras RS, Russell CM, Denys FR

Evaluation of the Syntac bonding system

Am J Dent 6 (1993) 17-21

Retief H, Mandras R, Russell CM

Shear bond strength required to prevent microleakage at the dentin/restoration interface

Am J Dent 7 (1994) 43-46

Rzanny A, Göbel R, Welker D

Quasistatische Prüfungen an 21 Füllungs- und Verblendkompositen

ZWR 104 (1995) 698-703

Rzanny A, Welker D

Biegefestigkeitsprüfungen an Kompositen im methodischen Vergleich

Quintessenz Zahntech 19 (1993) 1467-1473

Sakaguchi RL, Peters MC, Nelson SR, Douglas WH, Poort HW

Effects of polymerization contraction in composite restorations

J Dent 20 (1992) 178-182

Saunders WP

The marginal microleakage of resin-retained bridges in association with existing composite and amalgam restorations

Dent Mater 6 (1990) 20-23

Scheufele P

Suitability of different antagonists for pre-clinical in-vitro wear tests

J Dent Res 79 (2000) 279

Setcos JC

Heliomolar radiopaque als Amalgamersatz? Eine Fünf-Jahres-Studie von James C. Setcos, Manchester

Phillip J 12 (1995) 93-95

Shinkai K, Suzuki S, Katoh Y

Effect of an adhesive bonding system on wear resistance of resin composite restorations

Quintessence Int 28 (1997) 687-693

Shinkai K, Suzuki S, Leinfelder K, Katoh Y

How heat treatment and thermal cycling affect wear of composite resin inlays

J Am Dent Assoc 125 (1994) 1467-1472

Sjörger G, Hdlund SO, Jonsson C

A 3-year follow-up study of preformed beta-quartz glass-ceramic insert restorations

Quintessence Int 31 (2000) 25-31

Staehele HJ, Ackermann J

Experimentelle Untersuchungen über die Dichtigkeit von okklusalen Kompositfüllungen im Vergleich zu okklusalen Metallinlays

ZWR 100 (1991) 648-653

Steinmetz MJ, Pruhs RJ, Brooks JC, Dhuru VB, Post AC

Rechargeability of fluoride releasing pit and fissure sealants and restorative composites

Am J Dent 10 (1997) 36-40

Suliman AA, Boyer DB, Lakes RS

Cups movement in premolars resulting from composite polymerization shrinkage

Dent Mater 9 (1993) 6-10

Suliman AA, Boyer DB, Lakes RS

Interferometric measurements of cusp deformation of teeth restored with composites

J Dent Res 72 (1993) 1532-1536

Tani Y, Goto H, Ida K

Wear of posterior composite resins

Dent Mater J 6 (1987) 165-174

Taylor DF, Bayne SC, Leinfelder KF, Davis S, Koch GG

Pooling of long term clinical data for posterior composites

Am J Dent 7 (1994) 167-174

Tjan AHL, Chan CA

The polishability of posterior composites

J Prosthet Dent 61 (1989) 138-146

Turner CW, Meiers JC

Repair of an aged, contaminated indirect composite resin with a direct, visible-light-cured composite resin

Oper Dent 18 (1993) 187-194

Versluis A, Douglas WH, Sakaguchi RL

Thermal expansion coefficient of dental composites measured with strain gauges

Dent Mater 12 (1996) 290-294

Von Beetzen M, Li J, Nicander I, Sundstrom F

Microhardness and porosity of Class 2 light-cured composite restorations cured with a transparent cone attached to the light-curing wand

Oper Dent 18 (1993) 103-109

Wassel RW, McCabe JF, Walls AWG

Wear rates of regular and tempered composites

J Dent 25 (1997) 49-52

Wassel RW, McCabe JF, Walls AWG

A two-body friction wear test

J Dent Res 73 (1994) 1546-1553

Wassel RW, McCabe JF, Walls AWG

Subsurface deformation associated with hardness measurements of composites

Dent Mater 8 (1992) 218-223

Wassel RW, McCabe JF, Walls AW

Wear characteristics in a two-body wear test

Dent Mater 10 (1994) 269-274

Watts DC, Cash AJ

Analysis of optical transmission by 400-500 nm visible light into aesthetic dental biomaterials

J Dent 22 (1994) 112-117

Willems G, Lambrechts P, Braem M, Celis JP, Vanherle G

A classification of dental composites according to their morphological and mechanical characteristics

Dent Mater 8 (1992) 310-319

Willems G, Lambrechts P, Braem MJA, Vuylsteke-Wauters M, Vanherle G

The surface roughness of enamel-to enamel areas compared with the intrinsic roughness of dental composites

J Dent Res 70 (1991) 1299-1305

Xu X, Burgess Jo, Al-Owyyed, Weathersby

Artificial caries with fluoride-releasing and non-fluoride-releasing restorative materials

J Dent Res 79(2000) 284

Zeppenfeld W

Frei modellierte Brücken mit Heliomolar - Ein Fallbericht

Phillip J 12 (1995) 105-107

Publications on Heliomolar HB

Dunn JR, Munoz CA, Bernal G, Torres J, Wilson A
Clinical evaluation of Heliomolar HB for posterior restorations
J Dent Res 80 (2001) 91

Lowe E
Heavy duty tips and tricks for using Heliomolar HB
Dental Town Magazine 0 (2001) 14-15

Neme AL, Pink FE, Aksu MN, Maxson BB, Kashani N
In-vitro microleakage in class II packable resin composite restorations
J Dent Res 80 (2001) 105

Heintze SD, Cavalleri A, Vogel K, Sele J
Microleakage of composite class II restorations in vitro
J Dent Res 80 (2001) 589

Jackson RD
Replacing failed amalgam fillings using Heliomolar HB direct adhesive resin
Dental Products Report 3 (2001) 70

This documentation contains an overview of internal and external scientific data (information). The documentation and the corresponding information have been prepared exclusively for in-house use and for the information of external partners of Ivoclar Vivadent. They are not intended for any other purpose. Although we assume that the information complies with the latest standard of technology, we did not check all of them and may thus not assume any responsibility for their accuracy, truthfulness, or reliability. Liability cannot be assumed for the use of this information, even if we obtain contrary information. The information is used at the sole risk of the reader. Information is made available 'as received' without explicit or implied warranty regarding suitability (without reservation) for any specific purpose.

The information is made available free of charge. Ivoclar and its partners cannot be held accountable for any direct, indirect, immediate, or specific damage (including but not exclusively damage resulting from lost information, loss of use, or costs resulting from gathering comparable information), or for penal damages, which result from the use or failure to use this information, even if we or our representatives were informed about the possibility of such damage.

Ivoclar Vivadent AG
Research & Development
Scientific Service
Bendererstrasse 2
FL – 9494 Schaan
Liechtenstein

Contents: Dr Urs Lendenmann
Issued: November 2001
