Indications: Onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure, partial dentures.

Composition: Au 29.8, Pt 9.4, Pd 16.1, Ag 3.9, Cu 15.7, Zn 1.0, Ni 6.8, Ga 1.5, Fe 0.2, Cr 0.1, Other <1.0.

Instructions for Use: Istruzioni d'uso, Gebrauchsinformation, Instrucciones de uso, Instruções de uso, Brüksanvisning, Brugsanvisning, Käyttöohjeet, Produkttinformate.

Yellow, gold-based dental casting alloy, Type 4

Ivoclar Vivadent Worldwide: Australia, Brazil, Canada, China, Colombia, France, Germany, India, Japan, Mexico, Poland, Russia Federation, Sweden, Switzerland, Turkey, USA, UK.

ISO 13485 Quality Management System Certified Made in USA. www.ivoclarvivadent.com

EN INSTRUCTIONS FOR USE

MODELLATION Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration.

SPRUNG Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy).

BURN-OUT The suggested burnout temperature: 650-760C/1200-1400F High heat temperature investment: 650-760C/1200-1400F

MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1.

Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use.

CASTING Temperature: 965-1025C/1770-1880F Casting Temperature: 965-1025C/1770-1880F

METAL PREPARATION Carefully divest and clean the object with Al2O3, glass beads, or a pickling agent (such as Prevox®).

HEAT TREATMENT Annealing: 620C/1150F for 10 minutes; quench immediately (water) Hardening: 330C/625F for 30 minutes; air cool.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/112F.

POLISHING After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS Recommended for onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts, implant superstructure, partial dentures.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

France: Ivoclar Vivadent SAS, B.P. 118, 74410 Saint-Jorioz, Tel. +33 450 88 64 00, Fax +33 450 88 91 52.

SPRUNG METHOD

RECOMMENDATIONS: DIRECT: single crowns, inlays and onlays. INDIRECT: multiple units and multiple single crowns.

INSTRUKTIONER: 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.

2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.

SUGGESTIONS: 1. Thickness and design of wax pattern: follow your manufacturer's instructions.

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