Bluephase N[®] MC

EN Instructions for Use

Page 2

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Página 35

Appendix

Page 68

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Page 13

TR Kullanma Talimatı

Sayfa 46

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Página 24

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CE 0123

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist

Made in Austria

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

ivoclar
vivadent
clinical

Table of Contents


1. Product overview	3
1.1 List of parts	
1.2 Operation	
2. Safety	4
2.1 Intended use	
2.2 Indication	
2.3 Signs and symbols	
2.4 Safety notes	
2.5 Contraindication	
3. Start-up	6
4. Operation	7
5. Maintenance and Cleaning	9
6. What if ... ?	10
7. Warranty / Procedure in case of repair	11
8. Product specifications	12

1. Product overview

1.1 List of parts

1.2. Operation

The light is switched on and off by operating the start button.

2. Safety

2.1. Intended use

Bluephase N MC is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended use also includes the observation of the notes and regulations in these Instructions for Use.

2.2. Indication

Bluephase N MC is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 430–490 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporaries, as well as luting materials for brackets and indirect restorations, such as ceramic inlays.

2.3. Signs and symbols

Contraindication

Symbols on the curing light

Double insulation
(device complies with safety class II)

Protection against electrical shock
(BF type apparatus)

Observe Instructions for Use

Observe Instructions for Use

Caution

The appliance must not be disposed of as urban waste. Information regarding the disposal of appliance can be found on the respective national Ivoclar Vivadent homepage.


Alternating current

Direct current

2.4. Safety notes

Bluephase N MC is an electronic device and a medical product which is subject to IEC 60601-1 (EN 60601-1) and EMC directives, as well as the 93/42/EEC Medical Device Directive. The appliance complies with the relevant EU regulations and is classified as an LED Class 2 product.

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes in these Instructions for Use have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

2.5. Contraindication

Materials whose polymerization range is activated outside the wavelength range of 430–490 nm. If you are not sure about certain materials, please ask the manufacturer of the corresponding product.

Use without light probe.

Do not use the appliance near flammable or combustible substances.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Usage and liability

- Bluephase N MC must only be employed for the intended use. Any other uses are contraindicated. Do not touch defective, open devices. Liability cannot be accepted for damage resulting from misuse or failure to observe the Instructions for Use.
- The user is responsible for testing Bluephase N MC for its use and suitability for the intended purpose. This is particularly important if other equipment is used in the immediate vicinity of the curing light at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent (see Accessories). The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light probe is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the handpiece during operation.

Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

Assumption of impaired safety

If it has to be assumed that safe operation is no longer possible, the device must be disconnected from the power supply to avoid accidental operation. This may be necessary, for example, if the apparatus is visibly damaged or no longer works correctly. The apparatus is only completely disconnected from the power supply if the power cord has been pulled out of the plug-in socket.

Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are generally sensitive to light, who take medicine against light sensitivity or take photosensitizing medication, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device and they should wear protective goggles (orange) that absorb light below 515 nm.

Heat development

As it is the case with all high-performance lights, the high light output results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible damage. Therefore, the stipulated polymerization time, particularly in areas near the pulp (adhesive 10 seconds) must be observed. Uninterrupted curing times of more than 40 seconds on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane, or skin should be strictly avoided. Polymerize indirect restorations at intermittent intervals of 20 seconds each or use external cooling with an air stream.

3. Start-up

Check the delivery for completeness and any possible transportation damage (see delivery form). If parts are damaged or missing, contact your local Ivoclar Vivadent partner.

Power pack

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply.

Connect the power cord with the power pack and then connect with the power supply.

Handpiece

Remove the handpiece from its packaging and clean the light probe and the handpiece (see Maintenance and Cleaning section). The light probe can be removed and reattached by slightly rotating it.

After that, mount the anti-glare protection on the light probe.

4. Operation

Disinfect contaminated surfaces of the curing light as well as light probes and anti-glare cones before each use. Furthermore, make sure that the light output emitted by the curing light permits adequate polymerization. For that purpose, inspect the light probe for contamination and damage and check the light intensity at regular intervals (e.g. with the Ivoclar Vivadent radiometer “Bluephase Meter”).

Selecting the curing time

Observe the Instructions for Use of the material used when selecting the curing time. The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use of the relevant material, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 9 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

Curing times		800 mW/cm ² ± 10%
Restorative materials	Composites • 2 mm ¹⁾ IPS Empress Direct / IPS Empress Direct Flow Tetric EvoCeram / Tetric EvoFlow Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
	• 4 mm ²⁾ Tetric EvoCeram Bulk Fill / Tetric N-Ceram Bulk Fill Tetric Basic White	15 s 20 s
	Heliomolar / Heliomolar HB / Heliomolar Flow All conventional composites	20 s
	Compomers ³⁾ Compoglass F / Compoglass Flow	40 s
Indirect restorations / cementation materials	Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁵⁾ Dual Cement ⁵⁾ / Variolink II ⁵⁾	per mm ceramic: 10 s per surface
Adhesives	AdheSE / ExcITE F / Tetric N-Bond / Tetric N-Bond Self-Etch / Te-Econom Bond	10 s
Temporary materials	Telio CS Link / Systemp.link Systemp.inlay/onlay / Telio CS Inlay/Onlay / Fermit / Fermit N Telio Stains Telio Add-On Flow	15 s per surface 10 s 15 s 20 s
Miscellaneous	Helioseal / Helioseal F / Helioseal Clear Heliosit Orthodontic Vivaglass Liner	20 s
	IPS Empress Direct Color / IPS Empress Opaque MultiCore Flow / Multicore HB	30 s 40 s

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case e.g. with dentin shades)

2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case e.g. with dentin shades)

3) Applies to a maximum layer thickness of 3 mm
4) Applies to light-curing (use of the base past only)
5) Applies to dual-curing

Start

The light is switched on and off by operating the start button. During curing a signal beep sounds every 10 seconds.

Acoustic signals

- Start (Stop) 1 beep
- After 10 seconds, 1 beep
- After 20 seconds, 2 beeps
- After 30 seconds, 3 beeps
- After 40 seconds, 4 beeps
- After 50 seconds, 5 beeps
- After 60 seconds, 1 beep, light switches automatically off.

If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again.

Light intensity

The light intensity is maintained at a consistent level during operation. The light intensity is calibrated to $800 \text{ mW/cm}^2 \pm 10\%$ in conjunction with the 10 mm light probe supplied with the curing light.

If another light probe than the one supplied is used, the indicated light intensity may be directly affected.

In a light probe with parallel walls (10 mm), the diameter of the light entry and the light emission window is the same. When using focussing light probes (10>8 mm light probe, Pin-Point light probe 6>2 mm), the diameter of the light entry is larger than that of the light emission window. The incident blue light is thus bundled to a smaller area. As a result, the emitted light intensity is increased.

Pin-Point light probes are suitable for spot-on polymerization, e.g. to fix veneers before the removal of excess. For complete polymerization, the light probe must be changed.

Measuring the light intensity

The light intensity of the new LED polymerization device and the 10 mm light probe supplied with it can be checked with e.g. a Bluephase Meter.

If the measured value does not correspond with the expected light intensity, proceed as follows:

- Clean the light sensor of the radiometer, if dirty.
- Remove the light probe and clean the light emission window of the handpiece with a cotton bud dipped in alcohol.
- Clean the light probe, if dirty (see Maintenance and Cleaning).
- If the light probe is damaged, replace it with a new one.

If these measures do not improve the results, please contact your Ivoclar Vivadent partner.

5. Maintenance and Cleaning

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light probe. Disinfect contaminated surfaces of the device and anti-glare cones (e.g. FD 366/Dürr Dental, Incidin Liquid/Ecolab) and sterilize the light probe before each use if disposable protective sleeves are not used. Make sure that no liquids or other foreign substances enter the handpiece, the handpiece holder and particularly the power pack during cleaning (risk of electrical shock). Disconnect the handpiece from the power source when cleaning it.

Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean soiled plastic parts with a soap solution.

Light probe

Before cleaning and/or disinfecting the light probe, it should be pretreated. This applies to both automated and manual cleaning and disinfection.

Pretreatment

- Remove gross contamination immediately after use or within 2 hours. For this purpose, thoroughly rinse the light probe under running water (for at least 10 seconds). Alternatively, use a suitable aldehyde-free disinfecting solution to prevent the fixation of blood.
- To remove contamination manually, use a soft brush or soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula if necessary. Do not use sharp or pointed objects, as they may scratch the surface.

Cleaning and disinfection

For cleaning, immerse the light probe in a cleaning solution and ensure that it is sufficiently covered with liquid (ultrasound or careful brushing with a soft brush may support the effect). A neutral-enzymatic cleaning agent is recommended. When cleaning and disinfecting, please make sure that the agents used are free of:

- organic, mineral and oxidizing acids (the minimum admissible pH value is 5.5)
- alkaline solutions (the maximum admissible pH value is 8.5)
- oxidizing agents (e.g. hydrogen peroxide)

Afterwards, remove the light probe from the solution and thoroughly rinse it under running fresh water (for at least 10 seconds). Cleaning in a thermal disinfectant is an effective alternative.

Sterilization

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose. The sterilization time (exposure time at sterilization temperature) is 4 minutes at 134 °C (273 °F); pressure should be 2 bar (29 psi). Dry the sterilized light probe using either the special drying program of your steam autoclave or hot air. The light probe has been tested for up to 200 sterilization cycles.

After that, check the light probe for damage. Hold it against light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light probe with a new one.

Disposal

The curing light must not be disposed of as urban waste. Dispose unserviceable polymerization lights according to the corresponding legal requirements in your country.

6. What if....?

Indicator	Causes	Error rectification
<p>LED lights up in red, the device beeps</p> 	<p>Electronic component of the hand-piece may be defective</p>	<p>Disconnect the power cord from the device and insert it again. If the error reappears, please contact your Ivoclar Vivadent partner.</p>
<p>Device cannot be switched off</p>	<p>Electronic component is defective</p>	<p>Disconnect the power plug and remove it from the device.</p>

7. Warranty / Procedure in case of repair

The warranty for Bluephase N MC is 2 years from the date of purchase. Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases.

This is particularly valid for:

- damage resulting from improper handling, especially incorrectly stored batteries (see Technical Data: Transportation and storage conditions),
- damage to components resulting from wear under standard operating conditions (e.g. battery),
- damage resulting from external influences, e.g. blows, drop to the floor,
- damage resulting from incorrect set-up or installation,
- damage resulting from connecting the unit to a power supply, the voltage and frequency of which does not comply with the ones stated on the rating plate, and
- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the entire device including all its components must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent together with the purchase documents. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Centre. In case of a defect that cannot be rectified, please contact your dealer or your local Service Centre (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

8. Product specifications

Light source	LED
Wavelength range	430 – 490 nm
Light intensity	800 mW/cm ² ± 10 %
Operation	continuous operation
Light probe	Ø 10 mm, black, disinfectable and autoclavable
Anti-glare cones	autoclavable
Signal transmitter	Acoustic after 10 seconds and every time the start button / time selection button is pressed
Dimensions of the handpiece (without light probe)	L = 260 mm, W = 42 mm, H = 120 mm
Weight of the handpiece (incl. power pack and power cord)	420 g
Operating voltage handpiece	5 VDC with power pack
Power pack	Input: 100–240 VAC; 50-60 Hz; 0.28 – 0.14A Output: 5 VDC / 2 A Manufacturer: Friwo Type: FW7401M/05
Operating conditions	Temperature +10 °C to +35 °C (+50 °F to +95 °F) Relative humidity 35% to 75% Ambient pressure 700 hPa to 1060 hPa
Dimensions of the handpiece holder	D = 205 mm, W = 150 mm, H = 85 mm
Weight of the handpiece holder	220 g
Transportation and storage conditions	<ul style="list-style-type: none"> – Temperature –20 °C to +60 °C (-4 °F to +140 °F) – Relative humidity 10% to 75% – Ambient pressure 500 hPa to 1060 hPa – The curing light has to be stored in closed, roofed rooms. Protect the device from severe jarring.
Delivery form	<ul style="list-style-type: none"> 1 Handpiece holder 1 Handpiece 1 Power pack with power cord 1 Light probe 10 mm, black 3 Anti-glare cones 1 Instructions for Use

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Appendix

Guidance and manufacturer's declaration – electromagnetic emission

The Bluephase N MC is intended for use in the electromagnetic environment specified be-low. The customer or the user of the Bluephase N MC should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Bluephase N MC uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Bluephase N MC is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	N/A

Guidance and manufacturer's declaration – electromagnetic immunity

The Bluephase N MC is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase N MC should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line – line ± 2 kV line – earth	± 1 kV line – line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of the Bluephase N MC requires continued operation during power mains interruptions, it is recommended that the Bluephase N MC be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic Field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

The Bluephase N MC is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase N MC should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the Bluephase N MC, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Bluephase N MC is used exceeds the applicable RF compliance level above, the Bluephase N MC should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Bluephase N MC.

^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the "Bluephase N MC"

The Bluephase N MC is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Bluephase N MC can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bluephase N MC as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.4 \sqrt{P}$	80 MHz to 800 MHz $d = 0.4 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.13	0.13	0.22
1	0.40	0.40	0.70
10	1.3	1.3	2.2
100	4.0	4.0	7.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG

Bendererstrasse 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 9795 9599
Fax +61 3 9795 9645
www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723
Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400
Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd.
2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 3399
Fax +57 1 633 1663
www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118
F-74410 Saint-Jorioz
France
Tel. +33 4 50 88 64 00
Fax +33 4 50 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Germany
Tel. +49 7961 889 0
Fax +49 7961 6326
www.ivoclarvivadent.de

Wieland Dental + Technik GmbH & Co. KG

Schwenninger Strasse 13
D-75179 Pforzheim
Germany
Tel. +49 7231 3705 0
Fax +49 7231 3579 59
www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 22 2673 0302
Fax +91 22 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 6113555
Fax +39 051 6113565
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6903 3535
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37
Seochodong, Seocho-gu,
Seoul 137-855
Republic of Korea
Tel. +82 2 536 0714
Fax +82 2 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,
Piso 14, Col. Napoles
03810 México, D.F.
México
Tel. +52 55 5062 1000
Fax +52 55 5062 1029
www.ivoclarvivadent.com.mx

Ivoclar Vivadent BV

De Fruittuin 32
2132 NZ Hoofddorp
Netherlands
Tel. +31 23 529 3791
Fax +31 23 555 4504
www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 9999
Fax +64 9 914 9990
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawła II 78
00-175 Warszawa
Poland
Tel. +48 22 635 5496
Fax +48 22 635 5469
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/
office 10-06
115432 Moscow
Russia
Tel. +7 499 418 0300
Fax +7 499 418 0310
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204
P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 1 293 8345
Fax +966 1 293 8344
www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road
#02-01 San Centre
Singapore 169877
Tel. +65 6535 6775
Fax +65 6535 4991
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

C/ Ribera del Loira nº 46, 5ª planta
28042 Madrid
Spain
Tel. + 34 913 757 820
Fax + 34 913 757 838
www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14
S-169 56 Solna
Sweden
Tel. +46 8 514 939 30
Fax +46 8 514 939 40
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi
Sakayik Sokak
Nisantas' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 0802
Fax +90 212 343 0842
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE
United Kingdom
Tel. +44 116 284 7880
Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us