Description

Tetric® N-Flow Bulk Fill is a flowable, light-curing radiopaque composite for direct restorations in posterior teeth. As its opacity increases during polymerization, Tetric N-Flow Bulk Fill is also suitable for discoloured tooth structure. It is applied in increments of up to 4 mm as an initial layer in Class I and II restorations. Tetric N-Flow Bulk Fill cures with light in the wavelength range of 400–500 nm.

Shades

Tetric N-Flow Bulk Fill is available in three universal shades (^{IV}A, ^{IV}B and ^{IV}W).

Composition

The monomer matrix is composed of monomethacrylates and dimethacrylates (28 wt%). The fillers include barium glass, ytterbium trifluoride and copolymers (71 wt%). Additives, initiators, stabilizers and pigments are additional ingredients (<1.0 wt%). The total content of inorganic fillers is 68.2 wt% / 46.4 vol%. The particle size of the inorganic fillers ranges between 0.1 µm and 30 µm with a mean particle size of 5 µm.

Indication

- As initial layer / first increment in Class I and II composite restorations in permanent teeth
- Restorations in deciduous teeth

Contraindication

Placement of Tetric N-Flow Bulk Fill is contra-indicated

- if a dry working field cannot be established or the prescribed application technique cannot be applied.
- if the patient is known to be allergic to any of the ingredients in Tetric N-Flow Bulk Fill.
- Tetric N-Flow Bulk Fill must not be used as a luting composite.

Side effects

In rare cases, components of Tetric N-Flow Bulk Fill may lead to sensitization. The product must not be used in such cases.

Interactions

Materials containing eugenol/clove oil may inhibit the polymerization of composite materials. Consequently, the application of such materials together with Tetric N-Flow Bulk Fill must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application**1. Shade selection**

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide. The shade of the composite will correspond with that of the shade tab after the change in opacity has taken place, i.e. after polymerization.

2. Isolation

Appropriate isolation is required, preferably with a rubber dam (e.g. OptraDam® Plus).

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare any sharp internal edges or angles or additional undercuts in caries-free areas. The geometry of the cavity is predetermined by the dimensions of the caries lesion or the old filling. Slightly bevel or round out the enamel margins using finishing diamonds (grit-size 25–40 µm). Subsequently, rinse the cavity with water spray to remove all the residue and dry with water- and oil-free air.

4. Pulp protection / Base

In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide-based preparation (e.g. ApexCal®) and subsequently covered with pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Application of the matrix / interdental wedge

Use a wrap-around matrix for cavities affecting the proximal area or a sectional matrix band and secure it with wedges.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the following bonding agent: Tetric® N-Bond Universal (using the self-etching protocol or in conjunction with phosphoric acid etching).

7. Application of Tetric N-Flow Bulk Fill

- Tetric N-Flow Bulk Fill can be applied in increments of up to 4 mm thickness.
- Tetric N-Flow Bulk Fill must be covered with a layer of methacrylate based universal or posterior composite (e.g. Tetric® N-Ceram/Tetric® N-Ceram Bulk Fill/IPS Empress® Direct). Processing and finishing should be conducted according to the respective Instructions for Use.
- Proximal contacts must be established by means of a matrix system. The filling material does not exert any pressure on the matrix band. The matrix band can be shaped with a suitable hand instrument (e.g. a balltype condenser) before and/or during light-curing.
- Complete polymerization requires a sufficiently long exposure time. For the recommendations regarding exposure time and light intensity see table 1.
- When using a metal matrix, additionally polymerize the composite material from the buccal or lingual/palatal aspect after removing the matrix, if no Bluephase N® polymerization light is used or the light probe cannot be ideally positioned, e.g. distant to composite or diverging scattering angle.

8. Finishing / Checking the occlusion / Polishing

After polymerization, remove excess material with suitable finishers or fine diamonds. Remove proximal excess with diamond finishers, tungsten carbide finishers, finishing strips or flexible finishing discs. Check the occlusion and articulation and grind in the necessary adjustments to prevent premature contacts or undesired occlusal paths on the surface of the restoration. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Tetric® N-Flow Bulk Fill

EN Instructions for Use

Light-curing, flowable resin-based dental restorative material

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Additional information

- In the case of repairs, additional Tetric N-Flow Bulk Fill can be directly applied to polymerized material. If the Tetric N-Flow Bulk Fill restoration has already been polished, the surface in need of repair should be roughened with diamond grinders and wetted with Heliobond before a new layer of Tetric N-Flow Bulk Fill can be applied.
- Tetric N-Flow Bulk Fill has to be applied at room temperature. Cold material may be difficult to dispense.
-  For single use only. If Tetric N-Flow Bulk Fill is applied directly in the mouth of the patient, the Cavifil or the application cannula of the syringe must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).
- Do not use oxidizing disinfectants to disinfect syringes and Cavifils.
- The use of any other application cannula may make the material hard to extrude.
- The recommended thickness of the increment is based on hardness profile measurements.

Warning

- Avoid contact of unpolymerized Tetric N-Flow Bulk Fill with skin/mucous membrane and eyes. Unpolymerized Tetric N-Flow Bulk Fill may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effects of methacrylates.

Shelf life and storage

- Storage temperature: 2–28 °C
- Close syringes/Cavifils immediately after usage. Exposure to light causes premature polymerization.
- Expiry date: see note on Cavifils, syringes and packages.
- Do not use Tetric N-Flow Bulk Fill after the expiry date.

Table 1

Unit Programm	Bluephase N	Bluephase NM	Bluephase NMC
High Power	10 s (1200 mW/cm ²)	15 s (800 mW/cm ²)	15 s (800 mW/cm ²)
Soft Start	15 s	---	---
Low Power	20 s (600 mW/cm ²)	---	---

Light intensity	Exposure time
≥ 500 mW/cm ²	20 s
≥ 1'000 mW/cm ²	10 s

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.