Instructions for Use

■ PRODUCT DESCRIPTION

Au-based dental casting / metal-ceramic alloy, Type 4

■ INDICATIONS*

Inlays, Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts / Cores, Bars, Attachments, Partial Dentures

■ WAXING/MODELLATION

Design the framework in full contour or in a reduced anatomic shape taking the planned veneer into consideration. For wax patterns use mechanical retention. Single crown minimum thickness 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the framework design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required measures to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

■ SPRING

Provide the modelled single-tooth restoration or bridge framework with sprues of a suitable size. In general the design of the connectors sprues, whether near shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned in the heat center. The connector sprues between the reservoir and the casting should be a maximum of 10 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

■ INVESTING

Use phosphate / gypsum-bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN-OUT

Recommended burn-out temperature: 800 °C

■ MELTING AND CASTING

Tools: Propan 0,15 bar, Oxygen 0,35 bar
Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean ceramic /ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended amount of used material is 1:1. Do not use flux.
Casting Temperature: 1250 °C

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■ OXIDATION

Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation procedure.
Temperature: 800 °C. Holding time: 5 min; Vacuum: No
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing.
Use the appropriate ceramic veneering material, following the manufacturer's instructions.
Highest recommended firing temperature: 950 °C

■ HEAT TREATMENT

Hardening: 550 °C for 15 min; bench cool

■ **SOLDERING AND LASER WELDING**
Soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly. Use flux sparingly.
Pre Solder: Porta Optimum Lot 940 Flux: High Fusing Bondal Flux
Post Solder / Flux: Porta Optimum Lot 710 Bondal Flux
Laser Welding Wire: Portadur IN

■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ **CONTRAINDICATIONS**
For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ **SIDE EFFECTS**
Individuals with known sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be aware of the possibility for dental alloy to affect MRI results and to disclose the presence of dental alloys to the MRI technologist prior to conducting a test.

■ **INTERACTIONS**
Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and /or suitable protective mask is advised!

■ STORAGE CONDITION

Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions. For safety reasons, this alloy cannot be accepted for images resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and for any purpose not explicitly stated in the Instructions. These regulations also apply if the material is used in conjunction with products of other manufacturers.

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■ OXIDATION

Das Gerüst auf dem Brennträger providing adequate support. To achieve a uniform result follow the oxidation procedure.
Temperature: 800 °C. Holding time: 5 min; Vacuum: No
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing.
Use the appropriate ceramic veneering material, following the manufacturer's instructions.
Highest recommended firing temperature: 950 °C

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