Instructions for Use Mode d'emploi Instrucciones de uso Instruções de Uso Kullanma Talimatı Инструкция

- Dual-curing / light-curing, resin-based dental luting material
- Composite dual / photopolymérisable pour le collage des restaurations dentaires
- Cemento de fijación con base en composite de polimeración dual y fotopolimerizable
- Compósito de cimentação fotopolimerizável e de polimerização dual
- Diş hekimliği için dual sertleşen / ışıkla sertleşen yapıştırma kompoziti
- Стоматологический фиксирующий материал светового /двойного отверждения

For dental use only.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

CE 0123

ivoclar
vivadent®
clinical

6-45763/0612/6-spr/REV0

Instructions for Use

Description

Variolink® N is a dual-/light-curing luting composite system for the adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations.

Variolink N may also be applied in the light-curing technique only (e.g. luting of veneers). For this purpose, either the merely light-curing Variolink N Base or the merely light-curing Variolink N LC is to be used. Thus, an appropriate version of Variolink N is available for the various clinical requirements.

Syntac® is a tried-and-tested dentin-enamel adhesive system for direct and indirect restorations. Syntac can be used in combination with both light- and self-curing composites. The light-curing bonding agent Heliobond must be used prior to the application of the composite.

ExcITE® F DSC is a dual-curing, filled, fluoride-releasing adhesive for enamel and dentin bonding in conjunction with the total-etch technique.

Shades

Variolink N is available in 4 base shades and 2 catalyst shades as well as 2 levels of viscosity:

	A-D
Base	
bleach XL	BL 1
transparent	T
white	A1
yellow (universal)	A3
Catalyst high/low viscosity:	
transparent	T
yellow (universal)	A3

We recommend the base with the low-viscosity catalyst paste for the placement of crowns and partial crowns, the

base in combination with the high-viscosity catalyst paste for inlays/onlays and veneers. The dentist, however, is free to chose the proper viscosity for each indication. Additionally, the purely light-curing Variolink N LC is available in 3 shades: Clear, +1 and +2.

Working time

About 3.5 min. at 37 °C/99 °F

Mixing ratio

Mix base paste and catalyst in a 1:1 ratio.

Composition

The monomer matrix of Variolink N is composed of bis-GMA, urethane dimethacrylate, and triethylene glycol dimethacrylate. The inorganic fillers are barium glass, ytterbium trifluoride, Ba-Al-fluorosilicate glass, and spheroid mixed oxide. Additional contents: initiators, stabilizers and pigments. The particle size is 0.04–3.0 µm. The mean particle size is 0.7 µm.

	Monomer % wt.	Filler % wt.	Filler % vol.
Base	26.3	73.4	46.7
Catalyst			
low viscosity	27.9	71.2	43.6
high viscosity	22.0	77.2	52.0

Variolink N LC is composed of dimethacrylates, silicon dioxide and ytterbium trifluoride. Additional contents are initiators, stabilizers and pigments.

- The total content of inorganic filler is approx. 40%vol.
- Filler particle sizes range from 0.04-3.0 µm

Syntac Primer is composed of triethyleneglycol methacrylate, polyethyleneglycol dimethacrylate, maleic acid and ketone in an aqueous solution.

Syntac adhesive is composed of polyethyleneglycol dimethacrylate and glutaraldehyde in an aqueous solution. Heliobond contains Bis-GMA, triethyleneglycol dimethacrylate, stabilizers and initiators.

ExcITE F DSC contains HEMA, dimethacrylate, phosphonic acid acrylate, highly dispersed silicone dioxide, initiators,

stabilizers and potassium fluoride in an alcohol solution. The ExcITE F DSC applicator is coated with initiators.

Indication

Adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations (inlays, onlays, partial crowns, crowns, bridges, veneers and endodontic posts).

Contraindication

The application of Variolink N is contraindicated:

- if a dry working field cannot be established, or if the stipulated technique cannot be applied;
- if the patient is known to be allergic to any of the ingredients of Variolink N, Syntac Primer/Adhesive, Heliobond or ExcITE F DSC.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to individual components have been reported.

Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, the application of materials containing these substances is to be avoided. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn can impair the curing process. Thus, the syringe must not be disinfected using oxidative agents. The syringe can be disinfected by wiping it with e.g. medical alcohol. Alkaline jet media (e.g. Airflow) impair the effect of the adhesive.

Application

For detailed information also refer to the Instructions for Use of the respective products.

1 Removal of the temporary restoration and cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or the prepared tooth with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyl® fluoride-free). Rinse with

water spray. Subsequently, dry with water- and oil-free air. Avoid over drying.

2 Try-in and adjustment of the shade of the restoration

To aid in achieving optimum final restoration results, the restoration can be tried in with Variolink N Try-In pastes at this stage. To avoid fracture of the restoration, do not check the occlusion at this stage. Subsequently, adjust contact points if necessary and polish ground surfaces with ceramic or composite polishers. After try-in, thoroughly remove the Try-In paste with water spray and dry the restoration with water- and oil-free air.

3 Isolation

When applying an adhesive luting technique with composites, absolute isolation of the operative area – preferably with a rubber dam such as OpraDam® Plus, alternatively with cotton rolls and saliva ejector – is required. If the dental hard tissue is contaminated with saliva, it needs to be cleaned again as indicated in section 1.

4 Pre-treatment of the restoration surfaces

4.1.1 Leucite reinforced glass-ceramic and feldspar ceramic restorations (e.g. IPS Empress® Esthetic, IPS Empress CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 seconds or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

4.1.2 Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 seconds or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

4.1.3 Composite or fibre-reinforced composite restorations

- Blast the restoration surfaces (according to the instructions of the manufacturer of the restorative material).

- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse with water spray and dry with oil-free air

4.2 Apply a primer agent (e.g. Monobond® N or Monobond-S) with a brush or microbrush to the pre-treated surfaces. Allow the material to react for 60 seconds. Subsequently, disperse with a strong stream of air.

Important! Contamination with saliva, blood or water must absolutely be avoided during the entire procedure. In case of contamination, the surfaces of glass-ceramic, lithium disilicate glass-ceramic and composite restorations should be cleaned with the universal cleaning paste Ivoclean. Subsequently, the primer agent (e.g. Monobond N or Monobond-S) needs to be applied again.

5 Pre-treatment of the cavity

5.1 Using a matrix band in the proximal regions inhibits accidental etching of the adjacent teeth and facilitates the subsequent removal of excess Variolink N. Hold matrix band in place with a suitable interdental wedge.

5.2 Apply N-Etch (37% phosphoric acid gel) to the prepared enamel and then flow the etchant onto the prepared dentin. The acid should be left on the enamel for 15–30 seconds and on the dentin for 10–15 seconds.

IMPORTANT! Do not etch dentin longer than 15 seconds. Dentin etching is optional and not mandatory if Syntac is used.

5.3 Remove all etchant gel with a vigorous water spray for at least 5 seconds. Excess moisture should be removed until a slightly wet and shimmering dentin surface remains (wet bonding). This can be done with either an air gun, a dry brush, a foam pellet or other lint-free absorbents.

IMPORTANT! Do not overdry the dentin!

6 Application of the adhesive (e.g. Syntac or Excite F DSC)

6.1 Syntac

Attention! Do not mix Syntac Primer and Syntac Adhesive, but apply them successively!

6.1.1 Apply Syntac Primer to the cavity with a disposable brush and rub it gently into the surface. Syntac Primer requires a contact time of at least 15 sec on the dentin.

6.1.2 Disperse and thoroughly dry Primer excess with the air syringe. Do not rinse!

6.1.3 Apply Syntac Adhesive and allow it to react for 10 sec. Then thoroughly dry with blown air. Do not rinse!

6.1.4 Apply Heliobond and blow to a thin layer.

6.1.5 Heliobond is light-cured together with the Variolink N luting composite.

6.2 Excite F DSC

6.2.1 Apply Excite F DSC to the enamel and dentin and agitate the adhesive on the prepared surfaces for at least 10 seconds. Make sure that all the cavity walls are completely covered.

6.2.2 Disperse Excite F DSC to a thin layer with a weak stream of air, thereby removing any excess.

IMPORTANT! Avoid pooling of the adhesive on the cavity floor, as this can compromise the accuracy of fit of the definitive restoration. A glossy appearance shows that the tooth surface is completely sealed.

6.2.3 Pre-curing of Excite F DSC is not necessary if dual-curing luting composites are used. If pre-curing of the adhesive is preferred, Excite F DSC should be polymerized for 10 seconds at a light intensity of more than 500 mW/cm² (e.g. Bluephase[®] C8, LOW mode, 800 mW/cm²).

7 Mixing of Variolink N Base and Catalyst

– Mix Variolink N Base and Catalyst in a 1:1 ratio on a mixing pad for 10 sec immediately before application.

– When mixing is too hard, air may be trapped in the paste. This may delay polymerization and compromise mechanical properties.

– Variolink N is a light-/dual-curing material and thus sensitive to blue light (operatory light, ambient light). Avoid intensive light during application.

– The working time of the mixed Variolink N is about 3.5 min at a temperature of 37 °C/99 °F.

– **Attention!** Variolink N LC is a purely light-curing component and cannot be used in the dual-curing mode/mixed with the catalyst.

8 Placement of the restoration with Variolink N

8.1 Inlays, onlays, crowns and bridges

– Apply the mixed Variolink N to the area of the restoration which is to be bonded.

– First, place restoration in situ with slight pressure and remove excess with a brush or another suitable instrument.

– When working on large restorations, make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).

– Increase pressure and maintain for several seconds. Remove the Variolink N excess thus forced out with a brush.

8.2 Veneers

– Apply the Variolink N LC or the Variolink N Base with a brush or a spatula and/or with the help of the application tip directly to the preparation and/or to the inner aspect of the restorations, if required.

– Seat the restoration exerting light pressure and remove gross excess with a suitable instrument (e.g. spatula, brush).

– Make sure to remove excess material within the processing time, especially in areas that are difficult to reach (proximal areas and gingival margins).

– Increase pressure and hold for a few seconds.

– Maintain the pressure and light-cure the central area of the veneer for 3-4 seconds with an output of approx. 650 mW/cm² (e.g. Bluephase C8, LOW mode) to tack the restoration into place.

– **IMPORTANT:** Do not light-cure any proximal or marginal regions directly!

– Remove further excess with a suitable instrument.

When doing this, take care not to wipe any material out of the marginal areas.

9 Polymerization

Polymerize Variolink N segment by segment. Begin with the proximal margins.

- If a curing light with a light intensity of at least 800 mW/cm² is used (e.g. bluephase C8, HIP mode, 800 mW/cm²), polymerize for 10 sec per mm of ceramic and segment.
- If a curing light with a light intensity of 500 – 800 mW/cm² is used (e.g. LEDition, 600 mW/cm²), polymerize for 20 sec per mm of ceramic and segment.

Notes:

As all composites, Variolink N is subject to oxygen inhibition. This means that the surface layer (approximately 100 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

10 Finishing of the restoration

- Apply finishing and polishing strips in the proximal regions – Check occlusion and function and make corrections if necessary.
- Polish restoration with polishers (Astropol®, OptraPol® or OptraFine) or discs.

Warning

Unpolymerized Variolink N pastes and Heliobond may cause slight irritations.

Syntac Primer is an irritant and highly flammable.

Syntac Adhesive contains glutaraldehyde and is detrimental to health if inhaled or swallowed.

Excite F DSC is an irritant. Avoid contact with the skin, mucous membrane or eyes. In case of contact with the eyes, immediately rinse with water and seek medical advice. If the material comes into contact with the skin,

rinse with copious amounts of water. In rare cases, skin contact can cause sensitization to methacrylates and glutaraldehyde. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Variolink N, Syntac Primer, Syntac Adhesive, Heliobond or Excite F DSC after expiry date.
- Cool storage (2–8 °C/36–46 °F) is indicated for the Variolink N catalyst pastes
- Store Variolink N Base paste, Syntac Primer/Adhesive, Heliobond and Excite F DSC at 2–28 °C/36–82 °F.
- Variolink N syringes should be closed immediately after use. Exposure to light causes premature polymerization.
- Seal Syntac Primer/Adhesive and Heliobond bottles with the original cap immediately after use.
- Do not use Syntac Primer if the material shows clear signs of change of consistency.
- Shelf life: see expiry date on label

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions for Use or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions for Use.

The Material Safety Data Sheet is available online at www.ivoclarvivadent.com

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.