

# LEDition®

EASY. CONVENIENT. ECONOMICAL.

ERGONOMIC DESIGN

600 MW/CM<sup>2</sup>

MAINTENANCE-FREE LED AND  
CONTINUOUS OPERATION



**ivoclar**  
**vivadent**  
passion vision innovation



The **LEDition** polymerization unit is based on the latest LED technology, is easy to handle and always ready to use.



#### Easy.

Simply switch on and light cure... A universal program and one button only enable uncomplicated operation. In addition, the LEDition is very comfortable to hold due to its ergonomic design and low weight.

#### Convenient.

Thanks to the built-in fan, overheating of the unit is prevented – the prerequisite for a smooth continuous operation. The light performance is also reliable: normally, the LEDition light features a light intensity of 600 mW/cm<sup>2</sup> and is suitable for all the light-curing materials in the wavelength range of 430 – 490 nm.

#### Economical.

Another convenient advantage of LEDition: Other than with halogen lights, the light source does not have to be regularly replaced – this technological advantage not only saves costs but also time.

#### Technical data

Light source	3 watt LED
Wavelength range	430 – 490 nm
Max. light intensity	typically 600 mW/cm <sup>2</sup> (min. 500   max. 900 mW/cm <sup>2</sup> )
Light probe	10 mm, black, autoclavable
Cooling system	active
Acoustic signal	acoustical after 10, 20, 40, 60 etc. seconds
Weight of handpiece	200 g
Dimensions of handpiece	220 x 35 x 160 mm (L x W x H)
Weight of handpiece holder	245 g
Dimensions of handpiece holder	200 x 120 x 70 mm (L x W x H)
Power supply	power pack 100 – 240 VAC, 50/60 Hz
Warranty	1 year

#### Ordering information

Delivery form LEDition	Order No.
LEDition 100 -240 V	602108

#### Light probe

Light probe 10 mm, black	613700
Anti-glare cone	551756



Descriptions and data constitute no warranty of attributes.  
Printed in Liechtenstein

© Ivoclar Vivadent AG, Schaan/Liechtenstein  
603367/0209/e/BVD

Manufacturer and distributor  
**Ivoclar Vivadent AG**  
Benderstr. 2  
FL-9494 Schaan  
Principality of Liechtenstein  
Tel +423 / 235 35 35  
Fax +423 / 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

  
**ivoclar**  
**vivadent**<sup>®</sup>  
passion vision innovation